

Instructions for Use: ChannelCheck™ Convenience Pack

Brand Name of Product	ChannelCheck™ Convenience Pack
Generic Name of Product	3-in-1 Residual Soil Test Kit
Product Code Number(s)	UCC-222CP, UCC-222TS, ATS-2-CTP
Intended Use	To check internal channels for residual protein, blood and carbohydrates.
Range of Applications for Product	Internal channels that comes in contact with protein, blood and/or carbohydrate during clinical use.
Key Specifications of Product	Sensitivity of Reagent Pads: <ul style="list-style-type: none"> • Carbohydrate $\geq 25\mu\text{g/mL}$ • Protein $\geq 30\mu\text{g/mL}$ • Hemoglobin $\geq 0.25\mu\text{g/mL}$

Shipping & Storage	
Shipping Conditions & Requirements	Avoid direct sunlight.
Storage Conditions	<ul style="list-style-type: none"> • Box should be kept closed. • Keep in a cool dry place out of direct sunlight.
Packaging Contents	<ol style="list-style-type: none"> 1. 50 UCC-222 Test Strips (individually packaged) 2. 1 ATS Control Vial 3. 50 Ziplock Bags for collection 4. 50 Pre-filled 10 mL Water Syringes 5. 1 Interpretation Guide 6. 1 IFU
Shelf Life	Convenience Pack has a shelf life of 18 months.

Instructions for Using Product	
Description of Use (s)	
Preparation for Use	<ul style="list-style-type: none"> • The UCC-222CP is packaged as 5 boxes of 10 bags. These are labeled 1 of 5-box, 2 of 5, 3 of 5, 4 of 5, 5 of 5. Separating the boxes allows one box to be use at a time, and the other boxes to be stored until needed. • Box 1 states “This box contains the control test, control test must be ran first.” The control test is the first bag when box 1 is opened that says *Control Test* must be run first. Under the control test are 9 ea. HPC-222CP, individually wrapped ChannelChecks™, a syringe, and a 2x3 zip lock. Boxes 2-5 have 10 of the HPC-222CP. • Testing is to be conducted after cleaning and prior to disinfection/sterilization. • Control Test: The first step is to check the performance with the included vial of control soil. This is only done with one (1) of the individually packaged test strips and one (1) ATS control vial that is included in the Convenience Pack. To test, remove the vial of dehydrated test soil from the package. The test vial holds enough lyophilized test soil to create a single milliliter of test soil. <ol style="list-style-type: none"> 1. Rehydrate Soil: To rehydrate, unscrew the cap from the vial, then add exactly 1 mL of water to the vial. Screw the cap back on the vial ensuring you have a tight seal. 2. Shake Vigorously: Shake the vial vigorously for at least one minute. Check the vial to make sure the soil has been completely rehydrated. 3. Retrieve a Single Test Strip: Retrieve a single ChannelCheck™ test strip from the pack. 4. Dip Test Strip into Vial: Dip the test into the vial for 5 seconds, making sure to completely immerse all three test pads into the solution. 5. Dab Side of Test Strip on Absorbent Pad: After 5 seconds, remove the test strip and dab the side of the moistened test pad on a clean, dry absorbent pad to wick off excess water. 6. Wait 5 Minutes: The reagents in the test pads require time to interact with the residual soil, so wait a complete 5 minutes before reading the results. 7. Compare Results to Control Color Chart: After 5 minutes, compare the results to the Control Color Chart. The colors of each test pad should closely

approximate the colors on the Control Color Chart found on the Interpretation Guide.

8. **Record Results:** On a log sheet, record the results of each pad.

Diagrams (drawings, pictures)

ChannelCheck™ Convenience Pack



Figure 1



Figure 2



Figure 3



Figure 4



Figure 5



Figure 6




Figure 7

Steps for Use of Product

1. Uncap Pre-filled 10 mL Water Syringe.
2. **Flush the Water Through Channel:** This is done by flushing the channel(s) of the item with the 10 mL pre-filled syringe of water followed by the 10 mL of air to facilitate flushing. Place the slip tip at the channel to be tested and use the plunger-rod to deliver the water to sample the channel. Refill with air to finish the sampling procedure. (**Figure 4**)
3. **Recapture Water in the Ziplock Bag:** Recapture the water in a clean Ziplock Bag (see Ziplock Bag collection instructions below in the **Additional Information** section). (**Figure 5**)

	<ol style="list-style-type: none"> 4. Dip Test Strip into Water: Dip the test strip into the recaptured water, ensuring that all three pads are completely immersed. Keep test strip immersed for 5 seconds. (Figure 6) 5. Dab Side of Test Strip: Then remove it from the water. Dab the side of the test strip on a clean, absorbent surface to wick away excess water. 6. Wait 90 Seconds: The reagents in the test pads require time to interact with the residual soil, so wait a complete 90 seconds before reading the results. (Figure 7)
Interpretation of Results	<ol style="list-style-type: none"> 1. Compare to Color Chart: Compare test strip to the No Residue Color Chart found on the Interpretation Guide to interpret results. 2. Interpret Results: If the colors on any pad deviate from the No Residue Color Chart, this indicates a soiled item and it should be re-cleaned and re-tested until test results match the No Residues Color Chart.
Contraindications of Test Results	<ul style="list-style-type: none"> • Residual peracetic acid based disinfectants may interfere with the carbohydrate and blood pads of the ChannelCheck™. • Oxidizing agents like chlorine or hypochlorite may produce a color change on the hemoglobin pad. These residues should not remain after an item has been properly rinsed. • Excess residual Intercept® (brand of Cantel Medical) detergent can cause a color change on the protein pad. It is advised that you ensure adequate rinsing to remove any excess detergent prior to testing with ChannelCheck™. • Sterile water with adhered foil lids should not be used because of the potential for a reaction with the carbohydrate pad.
Documentation	Record Results: On a log sheet, record the results of each pad.
Special Warnings and Cautions	<ul style="list-style-type: none"> • Be sure to use the color chart that comes with the ChannelCheck that is included with this package. • ChannelCheck™ does not ensure that an item is safe for use, or free of contamination. It is to be one step in a total quality process implemented by the healthcare facility to verify the cleaning process. • Do not swirl the strip when dipping in the water. Swirling can cause color to run off the pad and change the results. • It is IMPORTANT that the test strips be protected from ambient moisture, light and heat to guard against altered reagent activity and deterioration. • It is possible some of the reagents in any one of the pads may be released when immersed in water, slightly coloring the water. This is normal and will not adversely affect the performance of the test.
Disposal	It is recommended to dispose of the used test strip in a suitable biohazard container.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	<p>Ziplock Bag Sample Collection:</p> <ol style="list-style-type: none"> 1. Open the plastic bag by gently pushing from the side of the bag. (Figure 1) (This will help create a wide enough opening so the clean plastic bag can be placed over the distal tip.) 2. Push the distal tip halfway down into the clean plastic bag. (Figure 2) 3. Once the tip is halfway into the clean plastic bag, seal the bag by pushing the sides together. Close the seal about ¾ of the way (up to the distal tip) and then stop. This will provide enough of a seal to capture the sample without the bag falling off during the sampling process. (Figure 3) 4. Follow the steps in Steps for Use of Product.

Related Healthmark Products	ATS2015
Other Product Support Documents	Cleaning Verification Brochure, Cleaning Verification Price List, ChannelCheck™ Specification Sheet, ChannelCheck™ Bottle Label, ChannelCheck™ Packaging Insert, ChannelCheck™ Validation Study, Instructions for Residual Soil Tests, Sample Policy with competency for ChannelCheck™, MSDS ChannelCheck™ UCC-222.
Reference Documents	<ul style="list-style-type: none"> • ALFA MJ, DEGAGNE P, AND OLSON N. WORST-CASE SOILING LEVELS FOR PATIENT-USED FLEXIBLE ENDOSCOPES BEFORE AND AFTER CLEANING. AM J INFECT CONTROL, 27:392–401, 1999. • ALFA MJ, DEGAGNE P, AND OLSON N. VALIDATION OF ATS AS AN APPROPRIATE TEST SOIL. ZENTR STERIL, 13(6):387–402, 2005. • ALFA MJ, OLSON N, DEGAGNE P, AND JACKSON M. A SURVEY OF REPROCESSING METHODS, RESIDUAL VIABLE BIOBURDEN AND SOIL LEVELS IN PATIENT-READY ENDO-SCOPIC RETROGRADE CHOLANGIOPANCREATOGRAPHY DUODENOSCOPES USED IN CANADIAN CENTERS. INFECT CONTROL HOSP EPIDEMIOL, 23:198–206, 2002.
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2020-08-06 Suzanne Latta