

# 7 ProChek™-W

Reorder# PT-SQ-001

- 1 Wear clean gloves. Open the Ziploc™ bag by gently pushing from the side of the bag (this will help create a wide enough opening so the clean plastic bag can be placed over the distal tip) and push the distal tip halfway down into the clean plastic bag.
- 2 Once the tip is halfway into the clean plastic bag, seal the bag by pushing the sides together. Close the seal about ¾ of the way (up to the distal tip) and then stop.
- 3 Flush the lumen with 10 mL of water. Collect the water in the provided Ziploc™ bag.
- 5 After recapturing the sample in the clean plastic bag, gently pull the clean plastic bag off the distal tip.
- 6 Unscrew the cap to the ProChek™-W vial. Pipette 1mL of sample from the ziploc™ bag.
- 7 Add 1 mL of sample to a ProChek™-W vial.
- 8 Replace the cap on the ProChek™-W vial.
- 9 Mix the contents in the vial by shaking at least 5 times. Wait 5 minutes.
- 10 Compare color of the solution to the provided interpretation chart.

## Instruction Guide for Residual Soil Tests




### How to Use this Booklet

- 1 Determine the test or test kit you are using.
- 2 Read and follow the instructions for the test to be used.
- 3 Record results on the log sheet.
- 4 Immediately report any test failure to department management.

RST-GUIDE: 03-2021

# 1 HemoCheck™

Reorder# HC-101


- 1 Open the pouch of the test kit. Included: A: indicator-vial (clear cap), B: activator-vial (green cap) and C: cotton swab. 
- 2 Wet surfaces are swabbed with the dry swab. Dry surfaces are swabbed by moistening the swab with a drop of clean chlorine-free water.
- 3 Open the indicator vial (A-clear cap) and transfer the liquid into the activator-vial (B-green cap).
- 4 Place the sample-swab (C) into the vial (head down into the liquid) and shake at least 5 times.
- 5 Check the swab within 30 seconds for a color change to blue-green, indicating blood residues on the tested surface. In the presence of a large amount of blood the whole indicator solution will change to dark blue. Record the result.



Positive for Blood

# 2 ProCheck™-II

Reorder# PT-313

- 1 Gloves must be worn to avoid test contamination. Included is one pouch with a purple capped vial with the test chemistry and a separately non-premoistened swab. Open pouch of the test kit. Moisten swab with the 3mL water vial that is provided. 
- 2 Swab the suspected target area.
- 3 Place swab into reagent vial. Recap the vial and shake at least 5 times.
- 4 If protein is present, color change of the liquid and/or on the swab to blue-green will occur. In case of soluble proteins, there will be an immediate color change. In the case of denatured proteins, color change can take up to 5 minutes.
- 5 Check liquid and swab for a color change to blue-green within 5 minutes. If no color change within five minutes the test is negative for protein. Record the result.



Blue-Green Positive for Protein

# 3

## HydroCheck™ Reorder# HYD-xxx

- 1 Insert the swab into the channel to be tested and feed it all the way through the channel until it exits the distal tip.
- 2 Immediately after the swab exits the distal tip, cut the swab approximately 2 cm behind the tip with scissors directly into the test vial.
- 3 Close and gently shake at least 3 times.
- 4 Fully reprocess the sampled item according to the manufacturer's IFU after testing with the HydroCheck™.
- 5 A color change to purple on the swab indicates residual moisture at or above the level of detection.
- 6 A small amount of residual moisture on the tip of the swab will be indicated by a small purple spot. Large amounts of moisture can turn the whole test purple.



# 4

## ProformanceQA.com

ProformanceQA™ is a secure online database that allows you to record key statistics for the performance of your decontamination procedures, including the test results from Healthmark's ProFormance™ line of monitoring products. With this cloud-based database you can record data, use the information to generate impactful reports and track the performance of your cleaning equipment over time.

To begin the setup process go to [www.proformanceqa.com](http://www.proformanceqa.com). Then, please call 800-521-6224 ext. 6657 or email [proformanceqa@hmark.com](mailto:proformanceqa@hmark.com) to schedule a call with the ProFormanceQA Help Desk.




# 5

## ChannelCheck™ Reorder# UCC-222

### Directions:

- 1 Using at least a 10ml syringe, fill with 10ml of prepackaged water and 10ml of air.
- 2 Flush lumen with pre-packaged water followed by air.
- 3 Recapture water in the provided zip-lock bag.
- 4 Dip test strip into water being sure that all three pads are completely submerged. Keep the strip submerged in the water for 5 seconds.
- 5 Remove strip from water and dab the side of strip on a clean, absorbent surface to wick away excess water.
- 6 Wait 90 seconds.
- 7 Compare test strip to the "No Residues" color chart to interpret results.
- 8 Should the color on any pad deviate from the "No Residues" this indicates a dirty instrument and it should be re-cleaned and re-tested.



# 6

## NOW! Test™ Reorder# NOW-1000

- 1 Don gloves and pick an item to test.
- 2 Flush lumen with pre-packaged sterile water, then purge with 30mL of air with an sterile syringe.
- 3 Recapture water in the provided zip-lock bag.
- 4 With the provided sterile pipette, draw up 0.5ml of water.
- 5 Add the 0.5 of water from previous step to the sterile cuvette with growth medium
- 6 Place vials in the block incubator and allow at least 12 hours of incubating.
- 7 Gently squeeze the bottom of the cuvette to help mix the reagent with the sample water.
- 8 Add 2 drops of Reagent A to the curvette and gently invert it a few times to help mix the reagent with the sample. Then turn on the fluorometer.
- 10 Remove the black cap from the fluorometer.
- 11 Place the cuvette in the fluorometer, line up the pointy side of the cuvette with the black line in the reader, and replace the black cap. Then on the screen, push "Measure."
- 13 On the screen, push "Blank."
- 14 On the screen, push "Measure." Then wait 10 minutes.
- 15 Any number above "200" indicates contamination and requires recleaning.