



# ProChek-II / PT-202 Instructions for use

## Protein Residue Instrument Assay:

The test kit for detection of protein residue on surfaces.



1. Gloves must be worn throughout the test procedure to avoid contamination of the test. Surfaces are swabbed with the supplied protein-free premoistened swab.
2. Swab the sample thoroughly. Concentrate swabbing action on difficult to clean areas like joints and crevices (visible spots or discolorations).
3. Unscrew cap of the vial. Do not touch the swab. Place the swab into the vial. Cut the swab end off into the vial with the supplied scissors. Place cap back on the vial and tighten and shake at least 5 times. If protein is present, color change of the liquid and/or on the swab to blue-green will occur. In case of soluble proteins, there will be an immediate color change. In the case of denatured proteins (often the case with instruments subjected to reprocessing) color change can take up to 5 minutes.
4. Check the liquid and the swab for a color change to blue-green within 5 minutes and no longer. If no color changes within 5 minutes the test is negative for protein. Record the result for quality assurance. 1 µg of denatured protein residue on the swab will develop a small blue-green spot. In the presence of large amounts of soluble protein the whole indicator solution will change blue-green.

NOTE: Positive controls can be used for verification and training purposes. If other swabs not supplied with the ProChek II-Test are used, a negative and positive control needs to be performed to exclude interferences, and protein free swab and water should be used.

### Marketing & Sales:

**Manufacturer:**  
PEREG GMBH  
Porschestr. 12  
D-84478 Waldkraiburg  
Germany

## -INFORMATION-

Make it visible!  
PEREG clean monitoring- the new perspective in CSSD

### PRINCIPLE

The formation of a protein-dye complex is used to detect small protein residue by means of a colour change on a swab. Swabbing is the preferred sampling technique in order to also detect insoluble residue.

### RANGE OF APPLICATION

Detection of protein residues on surfaces. Examples: Chamber walls of WD's, Ultrasonic cleaner, work bench, surfaces of surgical instruments or inside channels of endoscopes. Test for residues of protein based test soils.

### MEASURING RANGE

The test kit can detect down to 1 µg of insoluble protein by showing a blue-green spot.

### INTERFERENCES

Contact to alkaline substances (larger amount of alkaline detergent) can trigger colour change. Quaternary ammonium salts (used in some disinfectants) will give a false positive result. Contact to bare hands can transfer protein particles onto the swab and may give a false positive result.

### CONTENTS OF PACKAGE

15 X single use tests for detection of protein residue on surfaces and inside hollow instruments.

**STORAGE:** Store ProChek-II in a cool place at 2°C- 25°C. Keep away from light and heat and do not freeze.

**SHELF LIFE:** See imprint

Updated SMK 1/20/2012 Rev B