



## Instructions for Use: ProChek™-W Semi-Quantitative Protein Test for Robotic Arms

<b>Brand Name of Product</b>	ProChek™-W Protein Test for Robotic Arms
<b>Generic Name of Product</b>	Protein Test for Robotic Arms
<b>Product Code Number(s)</b>	PT-SQ-001
<b>Intended Use</b>	Testing for residual protein remaining on robotic arms after cleaning.
<b>Range of Applications for Product</b>	Robotic arms likely to be contaminated with protein soils during clinical use.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>● Measuring range for the test is 10 µg/ml-40 µg/ml or greater of protein</li> <li>● 15 Protein Reagent Vials</li> <li>● 15 Pipettes</li> <li>● 15 Zip lock Bags</li> <li>● 1 Interpretation Chart</li> </ul> <p>Customer will supply the syringe and prepackaged sterile water</p>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	Room temperature
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use(s)</b>	To test for residual protein that remain on robotic arms after cleaning.
<b>Preparation for Use</b>	<ul style="list-style-type: none"> <li>● Testing is typically conducted after cleaning prior to disinfection/sterilization.</li> <li>● Clean gloves</li> <li>● New 6 cc or larger syringe with blunt nose that fits the ports of robotic arms.</li> <li>● Customer will supply the syringe and prepackaged sterile water.</li> </ul>
<b>Diagrams of Use(s)</b>	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  <p><b>Figure 1</b></p> </div> <div style="text-align: center;">  <p><b>Figure 2</b></p> </div> </div>



**Figure 3**



**Figure 4**



**Figure 5**



**Figure 6**




**Figure 7**



**Figure 8**



**Figure 9**



**Figure 10**




**Figure 11**



**Figure 12**

1. Flush the lumen with 30cc of sterile wa
2. Collect the water in the provided Ziplo
3. Remove the can from the ProCheck™, 3

	 <p><b>Color of the solution when there is 10 µg, 20 µg, 30 µg, 40 µg of protein per 1 ML of the extract.</b>  <b>Note: For color accuracy please compare vials to printed interpretation guide included with the ProChek-W.</b></p>
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Wear clean gloves.</li> <li>2. Place zip lock bag on distal tip and secure it to ensure it does fall off.</li> <li>3. Fill syringe with 6 mL of sterile water. <b>Fig. 1</b></li> <li>4. With syringe, flush the lumen through Port 1. <b>Fig. 2</b></li> <li>5. Draw the water back out with the same syringe. <b>Fig. 3</b></li> <li>6. Gently eject the water into the supplied zip lock bag. <b>Fig. 4</b></li> <li>7. Make sure water comes in contact with entire area of distal tip and gently swish the distal tip in the recaptured water in the same zip lock bag for 10 seconds. <b>Fig. 5</b></li> <li>8. Remove the robotic arm from the zip lock bag. <b>Fig. 6</b></li> <li>9. Unscrew the cap to the ProChek™-W vial. <b>Fig. 7</b></li> <li>10. With supplied pipette, draw up 1 mL of recaptured water from the zip lock bag. <b>Fig. 8</b></li> <li>11. From the pipette, add 1 mL of water to a ProChek™-W vial. <b>Fig. 9</b></li> <li>12. Replace the cap on the ProChek™-W vial. <b>Fig. 10</b></li> <li>13. Mix the contents in the vial by shaking at least five times. <b>Fig. 11</b></li> <li>14. Wait 5 minutes.</li> <li>15. Place the ProChek vial in between the colors on the chart which will allow you to compare the color of the solution to the provided color interpretation chart. <b>Fig. 12</b></li> </ol> <ul style="list-style-type: none"> <li>● <b>The color interpretation chart shows the color of the solution when there is ≥ 10 µg, 20 µg, 30 µg, 40 µg or greater of protein per 1 ML of the extract.</b></li> </ul>
<b>Interpretation of Results</b>	<ul style="list-style-type: none"> <li>● In the case of soluble proteins, there will be an immediate color change.</li> <li>● In the case of denatured proteins (often the case with instruments subjected to reprocessing), color change can take up to 5 minutes.</li> <li>● For an accurate semi-quantitative result, wait the full 5 minutes.</li> </ul>
<b>Contraindications of Test Results</b>	Contact to alkaline substances (larger amount of alkaline detergent) can trigger color change. Quaternary ammonium salts (used in some disinfectants) will give a false positive result.
<b>Documentation</b>	Record test results on log sheet located at <a href="http://www.hmark.com">www.hmark.com</a>
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>● Evaluate the result immediately after 300 seconds (5 minutes) — late color changes are not valid.</li> <li>● Report a positive test result immediately.</li> <li>● A positive result is proof of remaining protein residue in the tested area</li> <li>● In case of a positive test result, the device should be reprocessed.</li> </ul>
<b>Disposal</b>	Since it is possible that organic soil is present, it is recommended to dispose all components of the ProChek™-W test kit in a biohazard container.

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	

<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	<ul style="list-style-type: none"> <li>• Lot control test results (Statement of Conformance) are available from Healthmark upon request.</li> <li>• Quality of Water for Testing: 1. It is recommended to use pre-packaged sterile water, such as sterile water for irrigation. 2. For effective testing, it is not necessary that the water remain sterile, but care should be taken not to contaminate after opening to avoid creating the opportunity for false positive test results.</li> </ul>
<b>Related Healthmark Products</b>	EndoCheck™, ProChek™, HemoCheck™, ChannelCheck™
<b>Other Product Support Documents</b>	ProFormance™ Brochure, ProFormance™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> hmark.com

2019-01-25 Suzanne Latta

Rev. B