


Instructions for Use: TempaChek™-LC2 Thermometer

Brand Name of Product	TempaChek™-LC2 Thermometer
Generic Name of Product	Liquid crystal thermometer
Product Code Number(s)	TEMP-LC2
Purpose of Product	To monitor the real time temperature of an instrument cleaning bath.
Range of Applications for Product	<ul style="list-style-type: none"> • Ultrasonic cleaners • Manual soaking baths • Sinks • Soak basins
Key Specifications of Product	<ul style="list-style-type: none"> • Provides real time feedback on the temperature of the bath • Adhesive allows placement of the thermometer below the waterline of the bath • Includes both Fahrenheit and Celsius scales: <ul style="list-style-type: none"> • 10-degree graduated scale from 60 to 150°F • 5- and 6-degree graduated scale from 16 to 65°C • Reusable • Precision to +/-1°C

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	1 per pack
Shelf Life	Expiration date printed on the package

Instructions for Using Product

Description of Use(s)	For monitoring the temperature of a cleaning bath in real time.
Preparation for Use	Be sure the surface for adherence is dry, clean and smooth.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Use one (1) TempaChek™-LC2 per bath. 2. Peel thermometer from release paper. 3. Apply to a clean area inside on the side wall of the tank below the intended fluid line and where the thermometer can be easily read. 4. TempaChek™-LC2 should be used to adjust the water temperature to the proper, optimal range for the cleaning agent used as specified by the manufacturer. 5. A color of green indicates the temperature is at the exact value. Blue is the temperature just below and tan is the temperature just above the actual temperature. Monitor the temperature during use of the bath to ensure compliance with target temperatures. Immediately report failure to achieve the target. 6. Record observations on a log sheet.
Interpretation of Results	Consult the detergent manufacturer's IFU to determine the proper temperature range for the cleaning bath. Temperature indicated by the TempaCheck™-LC2 should be in this range.

Contraindications of Test Results	
Documentation	Record observations on a log sheet.
Special Warnings and Cautions	
Disposal	<ul style="list-style-type: none"> • If adhesive residue remains on the surface, it may be removed with acetone or isopropyl alcohol solvent. • May be disposed of in regular trash.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	Proformance Brochure, ProFormance Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2020-08-06 Ralph J Basile