


## Instructions for Use: TempaChek™-LC Thermometer

|  |   |
|--|---|
| <b>Brand Name of Product</b>             | TempaChek™-LC Thermometer   |
| <b>Generic Name of Product</b>           | Liquid crystal thermometers   |
| <b>Product Code Number(s)</b>            | TEMPACHEK-LC, TEMP-LC2  |
| <b>Purpose of Product</b>                | To monitor the real time temperature of an instrument cleaning bath.  |
| <b>Range of Applications for Product</b> | <ul style="list-style-type: none"> <li>• Ultrasonic cleaners</li> <li>• Manual soaking baths</li> <li>• Sinks</li> <li>• Soak basins</li> </ul>   |
| <b>Key Specifications of Product</b>     | <ul style="list-style-type: none"> <li>• Liquid crystal thermometer provides real time feedback on the temperature of the bath.</li> <li>• Adhesive allows placement of thermometer below the water line of the bath.</li> <li>• 10 degree graduated scale from 60° to 150°F is easy to read. (Temp-LC2 include both Fahrenheit and Celsius scales)</li> <li>• Reusable.</li> <li>• Precision to +/-1°C.</li> </ul> |

### Shipping & Storage

|   |                                     |
|---|-------------------------------------|
| <b>Shipping Conditions &amp; Requirements</b> |                                     |
| <b>Storage Conditions</b>                     |                                     |
| <b>Packaging Conditions</b>                   | 1 per pack                          |
| <b>Shelf Life</b>                             | Expires one year from invoice date. |

### Instructions for Using Product

|  |  |
|--|--|
| <b>Description of Use (s)</b>            | For monitoring the temperature of a cleaning bath in real time.  |
| <b>Preparation for Use</b>               | Be sure surface for adherence is dry, clean and smooth.  |
| <b>Diagrams (drawings, pictures)</b>     |  <p style="text-align: center;"><i>TempaChek™-LC</i></p>   |
| <b>Steps for Use of Product</b>          | <ol style="list-style-type: none"> <li>1. Use one (1) TempaChek™-LC per bath.</li> <li>2. Peel thermometer from release paper.</li> <li>3. Apply to the inside, side wall of the tank below the intended fluid line and where the thermometer can be easily read.</li> <li>4. TempaChek™-LC should be used to adjust the water temperature to the proper, optimal range for the cleaning agent used as specified by the manufacturer.</li> <li>5. A color of green indicates the temperature is at the exact value. Blue is the temperature just below and tan is the temperature just above the actual temp.</li> <li>6. Monitor the temperature during use of bath to insure compliance with target temperatures. Immediately report failure to achieve the target.</li> <li>7. Record observations on log sheet.</li> </ol> |
| <b>Interpretation of Results</b>         | Consult the detergent manufacturer's IFU to determine the proper temperature range for the cleaning bath. Temperature indicated by the TempaCheck™-LC should be in this range.   |
| <b>Contraindications of Test Results</b> |  |
| <b>Documentation</b>                     | Record observations on log sheet.  |
| <b>Special Warnings and Cautions</b>     |  |
| <b>Disposal</b>                          | <ul style="list-style-type: none"> <li>• If adhesive residue remains on the surface, may be removed with acetone or isopropyl alcohol solvent.</li> <li>• May be disposed of in regular trash.</li> </ul>  |

### Reprocessing Instructions

|  |  |
|--|--|
| <b>Point of Use</b>                    |  |
| <b>Preparation for Decontamination</b> |  |
| <b>Disassembly Instructions</b>        |  |

|   |  |
|---|--|
| <b>Cleaning – Manual</b>                    |  |
| <b>Cleaning – Automated</b>                 |  |
| <b>Disinfection</b>                         |  |
| <b>Drying</b>                               |  |
| <b>Maintenance, Inspection, and Testing</b> |  |
| <b>Reassembly Instructions</b>              |  |
| <b>Packaging</b>                            |  |
| <b>Sterilization</b>                        |  |
| <b>Storage</b>                              |  |
| <b>Additional Information</b>               |  |

|  |  |
|--|--|
| <b>Related Healthmark Products</b>     |  |
| <b>Other Product Support Documents</b> | Proformance Brochure, ProFormance Price List   |
| <b>Reference Documents</b>             |  |
| <b>Customer Service Contact</b>        | Healthmark Industries Company, Inc<br>33671 Doreka<br>Fraser, MI 48026<br>1-586-774-7600<br><a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a><br><a href="http://hmark.com">hmark.com</a> |

2018-10-23 Ralph J Basile