


## Instructions for Use: TempaChek™-LC Thermometer

<b>Brand Name of Product</b>	TempaChek™-LC Thermometer
<b>Generic Name of Product</b>	Liquid crystal thermometers
<b>Product Code Number(s)</b>	TEMPACHEK-LC, TEMP-LC2
<b>Purpose of Product</b>	To monitor the real time temperature of an instrument cleaning bath.
<b>Range of Applications for Product</b>	<ul style="list-style-type: none"> <li>• Ultrasonic cleaners</li> <li>• Manual soaking baths</li> <li>• Sinks</li> <li>• Soak basins</li> </ul>
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Liquid crystal thermometer provides real time feedback on the temperature of the bath.</li> <li>• Adhesive allows placement of thermometer below the water line of the bath.</li> <li>• 10 degree graduated scale from 60° to 150°F is easy to read. (Temp-LC2 include both Fahrenheit and Celsius scales)</li> <li>• Reusable.</li> <li>• Precision to +/-1°C.</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	1 per pack
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use (s)</b>	For monitoring the temperature of a cleaning bath in real time.
<b>Preparation for Use</b>	Be sure surface for adherence is dry, clean and smooth.
<b>Diagrams (drawings, pictures):</b>	 <p style="text-align: center;"><i>TempaChek™-LC</i></p>
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Use one (1) TempaChek™-LC per bath.</li> <li>2. Peel thermometer from release paper.</li> <li>3. Apply to the inside, side wall of the tank below the intended fluid line and where the thermometer can be easily read.</li> <li>4. TempaChek™-LC should be used to adjust the water temperature to the proper, optimal range for the cleaning agent used as specified by the manufacturer.</li> <li>5. A color of green indicates the temperature is at the exact value. Blue is the temperature just below and tan is the temperature just above the actual temp.</li> <li>6. Monitor the temperature during use of bath to insure compliance with target temperatures. Immediately report failure to achieve the target.</li> <li>7. Record observations on log sheet.</li> </ol>
<b>Interpretation of Results</b>	Consult the detergent manufacturer's IFU to determine the proper temperature range for the cleaning bath. Temperature indicated by the TempaCheck™-LC should be in this range.
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	Record observations on log sheet.
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	<ul style="list-style-type: none"> <li>• If adhesive residue remains on the surface, may be removed with acetone or isopropyl alcohol solvent.</li> <li>• May be disposed of in regular trash.</li> </ul>

Reprocessing Instructions	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	

<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Proformance Brochure, ProFormance Price List
<b>Reference Documents</b>	
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2017-12-22 Ralph J Basile