

Instructions for Use: LTA MicroFlex Brush

Brand Name of Product	LTA MicroFlex Brush	
Generic Name of Product	Channel Brushes	
Product Code Number(s)	730-20, 730-50, 740-30, 740-50, 740-75, 750-30, 750-50, 750-75, 760-30, 760-50, 760-75, 770-30, 770-50, 770-75, 780-30, 780-50, 780-75	
Intended Use	For cleaning internal channels.	
Range of Applications for Product		
Key Specifications of Product	Components for item code 740 and 750	Material
	Flexible wire	Metallic Alloy
	Wire coating	Thermoplastic
	Connection	Stainless steel
	Protection cap	Stainless steel
	Fibers	Thermoplastic Compound with a low absorbency level
	Wire	Stainless steel
	Components for item code 760-770-780	Material
	Soul flexible stem	Metallic Alloy
	Spring	Thermoplastic
	Connection	Stainless steel
	Spring	Stainless steel
	Protection cap	Stainless steel
	Fibers	Thermoplastic Compound with a low absorbency level
Wire	Stainless steel	

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	
Shelf Life	

Instructions for Using Product	
Description of Use(s)	For the cleaning of internal channels.
Preparation for Use	<ul style="list-style-type: none"> It is important to match the diameter and length of the brush you are using with the channel you are cleaning. The Healthmark brush-sizing chart or similar type of tool may be used to select the proper diameter brush. Ensure the material composition of the brush bristles is suitable for the item to be cleaned. This should be in compliance with the item manufacturer's IFU. Brushes should be inspected prior to each use to be certain they are not worn, frayed or damaged. If discrepancies are noted, the brush should be discarded and not used. Brushes should also be inspected between uses for any residual organic soil that may cross contaminate the next item cleaned. If there is visual soil, the brush should be cleaned in accordance with the method(s) described below.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ul style="list-style-type: none"> Consult the manufacturer's IFU for the item you are cleaning and follow those instructions closely. Insert the brush by non-brush end into the working channel and then follow with the bristles. Push until the brush exits the other side of the instrument. This insertion technique optimizes the cleaning action and prevents/reduces kinking or bending at the point where the brush is inserted.

	<ul style="list-style-type: none"> • Then pull the brush in a traction motion. Brushes can be used for pushing when cleaning as well, however the traction method is recommended in order to get easy cleaning and avoid retro-contamination of the working channel. • It is important to make sure the brush is long enough to exit the lumen or channel completely when brushing so one can examine the bristles and clean them before used again. • While brushing, the brush should pass freely through the lumen. If significant resistance is encountered, stop advancing the brush and pull back - continuing to brush through the resistance might cause damage to the lumen.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> • Tube style brushes are intended for use as cleaning tools only. They are not for use in patient care or aseptic techniques. • When inserting in channel, do not hold brush at compromising angle that would cause it to bend. • Not to be used for cleaning plastic and/or fragile materials such as silicone. • Always protect yourself from brush splatter and contamination. • Clean between each use and disinfect at least daily, preferably between each use in order to limit the chance for cross contamination. See instructions for reprocessing below.
Disposal	Because MicroFlex brushes are used to clean patient used medical devices, it is recommended to dispose of them in the biohazard container with other contaminated devices.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	<ol style="list-style-type: none"> 1. Rinse gross contaminants from the brush in a deep sink and a spray apparatus, such as a shower gun. 2. Use the facility's cold water supply for this rinsing.
Disassembly Instructions	
Cleaning – Manual	<ul style="list-style-type: none"> • Brushes may be cleaned manually utilizing a high alkaline instrument detergent or an enzymatic detergent with the suitable enzymes for cleaning the types of soil the brush may come in contact with during cleaning (i.e., protease, lipase, etc.). • The IFU for the cleaning agent should be followed closely to ensure sufficient exposure time is allowed for the cleaning solution to work. • It is recommended that agitation accompany the cleaning process to improve results.
Cleaning – Automated	<ul style="list-style-type: none"> • Machine cleaning is recommended in a washer disinfector cleared by the FDA. • Products should be positioned in the washer to allow maximum water penetration and drainage. • No overlapping - partially covered devices will not be washed properly • Brushes should be secured in some way, such as the use of a wire or perforated tray that has a wire or perforated lid to prevent brushes from blowing around the wash chamber during cleaning. A container, such as the brush cleaning tube from Healthmark, is a suitable solution. • Products can be cleaned with alkaline, acidic and neutral detergents. • For final rinse, DI or RO water is recommended. • Use the “instrument” cycle of the washer-disinfector to clean the brushes
Disinfection	<ul style="list-style-type: none"> • The thermal disinfection stage of an automated washer is sufficient for disinfection of the cleaning brushes. • Brushes may be disinfected with liquid chemical sterilants in accordance with the disinfection manufactures' IFUs. Be sure to check material compatibility information from the disinfectant manufacturer. Should be compatible with nylon and with stainless steel.
Drying	
Maintenance, Inspection, and Testing	<ol style="list-style-type: none"> 1. Inspect for wear, fraying or damaged bristles or the twisted wire shaft. If observed, the brush should be discarded in compliance with the disposal instructions described above. 2. Brushes should also be inspected between uses for any residual organic soil that may cross contaminate the next device cleaned. 3. Maximum recommended use life: 3 months
Reassembly Instructions	

Packaging	
Sterilization	<ul style="list-style-type: none"> • Terminal sterilization of brushes is not normally necessary. However, these brushes are compatible with steam sterilization. • 740: maximum number of recommended sterilization cycles is 30 cycles. • 750: maximum number of recommended sterilization cycles is 60 cycles. • If steam sterilizing, recommended cycles are: 132°C/4 minutes or 135°C/3 minutes. • Maximum time/temperature sterilization cycle is 30 minutes at 135° C
Storage	Avoid compression during sterilization and storage
Additional Information	Follow all manufacturers' directions for proper usage and disposal of all cleaning agents and chemical disinfectants agents.

Related Healthmark Products	Brush Cleaning Tools
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2020-12-01 Ralph J Basile