

Instructions for Use: Sample Extraction Tools

Brand Name of Product	Sample Extraction Tool
Generic Name of Product	Single-use sample extraction tool
Product Code Number(s)	CKB-200-200, CKB-230-400, CKB-230-600
Intended Use	The sample extraction tool is intended for extracting samples to test for bacterial growth.
Range of Applications for Product	Extracting samples for testing
Key Specifications of Product	<ul style="list-style-type: none"> ● Proper length and diameters for contact with the internal surfaces. <ul style="list-style-type: none"> - CKB 200-200: 200cm Length x 2mm Diameter - CKB 230-400: 230cm Length x 4mm Diameter - CKB 230-600: 230cm Length x 6mm Diameter

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use(s)	To extract samples for testing for bacterial growth.
Preparation for Use	<ol style="list-style-type: none"> 1. Select the diameter and length for the internal surface/channel. 2. Remove the extraction tool from the packaging.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Insert the sample extraction tool into the channel. 2. Continue to feed the extraction tool through the length of the channel until it reaches the distal end. Never force the extraction tool, to avoid possible damage. 3. Once the extraction tool exits the distal end, cut the entire head of the bristle portion of the extraction tool into the sample collection container. 4. Carefully pull the entire sample extraction tool completely out of the distal end. 5. Dispose of the extraction tool after use.
Interpretation of Test Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> ● Select the extraction tool based upon the diameter and length of the channel chosen. ● Do not force the extraction tool into a channel to avoid possible damage.
Disposal	Treat used sample extraction tools as a biohazard and dispose of in accordance with facility policy of contaminated items.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	

Additional Information	
Related Healthmark Products	Flexible Endoscope Sampling Kits
Other Product Support Documents	ProFormance™ Brochure, ProFormance™ Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2019-01-30 Suzanne Latta

Rev. A