



<b>Brand Name of Product</b>	Wash Cradle
<b>Generic Name of Product</b>	Instrument wash cradle
<b>Product Code Number(s)</b>	IWC-001
<b>Purpose of Product</b>	The Wash Cradle is designed to hold surgical instruments in the open position to promote effective cleaning in the automated washer with a minimum amount of staff effort.
<b>Range of Applications for Product</b>	For hinged instruments that need to be held open during the washing process.
<b>Key specifications of Product</b>	<ul style="list-style-type: none"> <li>• 11.5” Length, 2” Width, 3.75” Height</li> <li>• Manufactured of stainless steel</li> <li>• Adjustable rack length</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	The design will increase effectiveness of the cleaning process of ring handled/hinged instruments by presenting them in an open and upright position. This positioning helps expose the most difficult to clean locations to washing (water & detergents) and instrument lubricants.		
<b>Preparation for Use</b>	Methods to Load: 1) Instruments are loose in basket.		
<b>Diagrams (drawings, pictures)</b>	<p><b>Assembly</b></p>		

**Loading**



7



8



9



10



11



12



13



14

**Steps for Use of Product**

**Assembly**

1. Cradle and key lock.
2. Install key lock by sliding open slit over wire guides.
3. Pinch wire guides together.
4. Turn the key lock until wire guides snap into place.
5. Slide key lock back.
6. Ready to load.

**Loading**

7. Be sure to open instruments to maximize exposure for cleaning.
8. Place each instrument onto the cradle rack.
9. Continue loading instruments onto the rack until all instruments are in place.
10. Slide key lock to place to secure instruments.
11. Pick up the wash cradle to place into basket.
12. Place the wash cradle into the basket.
13. Be sure that the wash cradle is secure in the basket. It can be used in various styles of baskets as long as the instruments are secure.
14. Instruments are ready to be washed.

**Interpretation of Results**

**Contraindications of Test Results**

**Documentation**

<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	Washer baskets, Nic Nac Washer Baskets & Balls
<b>Other Product Support Documents</b>	Prosys™ Brochure, Prosys™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2018-05-16 msmith