





**Instructions for Use: Weisenbach Instrument Forceps Stringer**

<b>Brand Name of Product</b>	Weisenbach Instrument Forceps Stringer
<b>Generic Name of Product</b>	Stringers
<b>Product Code Number(s)</b>	10-1775, 10-1776
<b>Intended Use</b>	Keeps ring handled instruments maximally exposed during sterilization.
<b>Range of Applications for Product</b>	Used for ring handled instruments for sterilization.
<b>Key Specifications of Product</b>	Manufactured of stainless steel.

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	N/A
<b>Storage Conditions</b>	N/A
<b>Packaging Contents</b>	N/A
<b>Shelf Life</b>	N/A

**Instructions for Using Product**

<b>Description of Use(s)</b>	Keeps ring handled instruments maximally exposed during sterilization.
<b>Preparation for Use</b>	<ol style="list-style-type: none"> <li>1. This requires two stringers to hold the rings open during sterilization.</li> <li>2. Each bar will hold one ring handle for that instrument.</li> </ol>

<b>Diagrams (drawings, pictures)</b>	 <p style="text-align: center;"><b>Figure 1</b></p>  <p style="text-align: center;"><b>Figure 2</b></p>
--------------------------------------	--



**Figure 3**

<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Unlock the “T” of the stringer by placing it in a vertical position. <b>Fig. 1</b></li> <li>2. Place the ring handled instruments in the open position.</li> <li>3. The rings of the instrument are put through both bars of the stringer. <b>Fig. 2</b></li> <li>4. When finished placing the instruments through the bars, the “T” is then locked back in the horizontal position. <b>Fig. 3</b></li> </ol>
<b>Interpretation of Test Results</b>	N/A
<b>Contraindications of Test Results</b>	N/A
<b>Documentation</b>	N/A
<b>Special Warnings and Cautions</b>	N/A
<b>Disposal</b>	N/A

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	N/A
<b>Preparation for Decontamination</b>	N/A
<b>Disassembly Instructions</b>	N/A
<b>Cleaning – Manual</b>	Before initial use of the stringer, cleaning is required. Manually wash stringer(s) with a neutral pH detergent and warm water.
<b>Cleaning – Automated</b>	Use the “instrument cycle” of the washer-disinfector compliant to ISO 15883-1 & 2.
<b>Disinfection</b>	The thermal disinfection stage of a washer-disinfector compliant to ISO 15883-1 & 2.
<b>Drying</b>	N/A
<b>Maintenance, Inspection, and Testing</b>	N/A
<b>Reassembly Instructions</b>	N/A
<b>Packaging</b>	N/A
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• See instructions under “<b>Steps for Use of Product</b>” for directions on how to string instruments.</li> <li>• When using for sterilization of instruments, open to the optimal width for use in the instrument set. Consult the instrument manufacturer’s IFU for specific instructions.</li> <li>• Stringers are materially compatible with steam and EO sterilization.</li> <li>• Steam sterilization at 132°C (270°F) for 4 minutes or 135°C (275°F) for 3 minutes.</li> <li>• Maximum steam sterilization at 135°C (275°F) for 30 minutes.</li> </ul>
<b>Storage</b>	N/A
<b>Additional Information</b>	N/A
<b>Related Healthmark Products</b>	Stringers
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	N/A
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyndale Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

