

## Use of the Flexible Inspection Scope (FIS) for Inspection of Flexible Endoscopes and Other Instruments

Product Code: FIS-005

Mary Ann Drosnock, MS, CIC, CFER, RM (NRCM), FAPIC



The Flexible Inspection Scope is a type of video borescope that includes a distal tip composed of a light source and camera lens at the end flexible shaft that comes in various lengths and sizes. A borescope is a piece of equipment used to inspect the inside of an instrument through a small opening (lumen/channel) found on that instrument. The Healthmark Flexible Inspection Scope is intended for inspection of any medical or surgical instrumentation preferably after cleaning, and prior to disinfection or sterilization. It may be used on the outside or inside surfaces of instruments including those such as arthroscopy routers, shavers, reamers, flexible endoscopes and other

lumened items.

According AAMI ST91 <sup>1</sup>, healthcare facilities should establish a comprehensive quality assurance and safety program to monitor all aspects of endoscope processing. This program should incorporate both visual inspections and testing of the equipment to identify conditions that may affect the cleaning or disinfecting processes.

Visual inspection alone may not be sufficient for assessing the efficacy of cleaning processes. Use of enhanced visualization methods, including inspection with a borescope, such as the Flexible Inspection Scope (FIS), is being incorporated into current standards and recommendations.

Visual inspection of the equipment should include the following:

- Residual organic soil
- Cracks and other damage to the instrument
- Integrity of fiber optic bundles
- Use of magnification

Both ST91 and AORN Endoscope Reprocessing Guidelines <sup>2</sup> state to consider inspection with borescope. According to ST91, “Inspection using magnification and additional illumination might identify residues more readily than the unaided eye... tools such as **video borescopes** of an appropriate dimension (length and diameter) may be used to visually inspect the internal channels of some medical devices...the use of methods that are able to quantitatively or chemically detect organic residues that are not detectable using visual inspection should be considered and included in facility policies and procedures on device cleaning...”

According to AORN VII.c.1., “Internal channels of flexible endoscopes may be inspected using an endoscopic camera or **borescope**. [2: High Evidence] Endoscopic cameras and borescopes penetrate the lumen and allow for improved visual inspection.”

In AAMI ST79 <sup>3</sup> Visual inspection is described as a verification of the Cleaning process. Section 7.6.4.5 states the following:

- After completing the cleaning process, personnel should visually inspect each item carefully to detect any visible soil.
- Inspection using magnification and light borescopes might identify residues not observable by the unaided eye.
- Visual inspection alone may not be sufficient for assessing the efficacy of cleaning processes; the use of methods that are able to measure organic residues that are not detectable using visual inspection should be considered in facility cleaning policy and procedures (see Annex D for available methods).


Also, the IAHCMM Endoscope Reprocessing Manual <sup>4</sup> from 2017 states the following: “The most important areas for inspect on are the lumens that run through the endoscope. Lumens pose a cleaning challenge because of their narrow structure that prevents visualization during the cleaning; ‘therefore, it is important to always check lumens for cleanliness after cleaning. Visual inspection of lumens can be accomplished using a borescope, a small flexible fiberoptic device that enables visualization of otherwise inaccessible areas within endoscope lumens.”

According to the SGNA POSITION STATEMENT on Management of Endoscopic Accessories, Valves, and Water and Irrigation Bottles in the Gastroenterology Setting<sup>5</sup>, “A comprehensive quality control program for reusable medical devices should be implemented and include, but not limited to the following:


- a. Visual inspection and equipment testing to identify conditions that may affect the cleaning or disinfection process. Damaged reusable items should be removed from use. Follow facility protocol for returning device.
- b. Procedures for monitoring the useful life of medical devices that include visual inspection, scheduled maintenance, and removal of equipment from use based on manufacturer's guideline.

**Areas to inspect:**

1. **Endoscopes:** There are many areas in an endoscope which can be inspected with the FIS including instrument suction channel, channel openings/valve housings, distal tip, forceps elevator recess on endoscopes. Here are some other suggested areas that might be considered for inspection on a flexible endoscope:



2. **Shavers:** Areas to inspect include the following:



## **IFU Support for using a Flexible Inspection Scope:**


1. Arthrex Hand Piece<sup>6</sup>: The IFU states the following, “Check the device for visible soil. It is recommended that the cannulation be inspected with an illuminated, magnifying scopes. Clean the device using the guidelines for manual cleaning if any soil is visible...”.
2. STRYKER Shaver Hand Piece<sup>7</sup>: The IFU states the following, “Visually inspect the hand piece, including all internal surfaces, for remaining soil. Use an endoscopic camera and endoscope if necessary to see the inner surface of the lumen. If soil remains, repeat the manual cleaning procedure, focusing on those areas...”.

In an FDA Safety Communication<sup>8</sup>, “Consider inspecting the inside of the devices following cleaning to ensure that they have been cleared of any tissue or fluids. There may be multiple ways to accomplish this. As one example, the facility that brought this situation to our attention uses a 3mm video scope to inspect the channels of the shaver handpiece.”

## **Study data:**

1. In a recent study conducted by Ofstead and Associates<sup>9</sup>, the Flexible Inspection Scope was used to inspect endoscopes in a healthcare facility. During this study there was good adherence to published reprocessing steps, however, upon inspection of the interior components of endoscopes with a FIS, borescopic examination revealed the following:
  - Residual fluid
  - Irregular surfaces and brown staining
  - Scratches, non-intact lining, and brown staining

Therefore, by examining inside the reprocessed endoscopes, the study revealed both damage and debris retained inside the channels of the inspected scopes. Inspections revealed damage and debris inside channels and residual fluid in channels and ports. These findings indicated that current reprocessing methods were not sufficient to result in properly reprocessed endoscopes. By identifying these findings, the facility was able to implement additional measures to ensure that endoscopes were being reprocessed to a higher standard of care.



2. In another publication by Ofstead and Associates <sup>11</sup>, the Flexible Inspection Scope was used to inspect GI endoscopes in a Healthcare facility over time. During the study, it was found that most GI endoscopes had irregularities upon inspection. Specifically, examination with the flexible inspection scope demonstrated scratches, discoloration, debris and fluid inside the endoscope channels.

Therefore, by performing inspection with this borescope, researchers were able to identify endoscopes that needed repair and had critical defects. These findings support use of enhanced visual inspection by use of a borescope, such as the Flexible Inspection Scope.

In a recent study entitled, “Inspection of Endoscope Instrument Channels After Reprocessing Using a Prototype Video Camera: A Pilot Study <sup>12</sup>,” authored by Thaker et al., it was determined that an inspection program with a video camera device showed that internal damage and/or discoloration of the instrument channel appears to occur frequently, even in newer endoscopes. Also, the study demonstrated that video inspection of the endoscope channel may be useful to audit reprocessing performance and to identify damaged endoscopes.

3. Ofstead and associates presented the study below at the AORN 2017 Conference<sup>12</sup>. In this study, it was found that sterilized ureteroscopes had high levels of visible

irregularities, debris, damage and contamination upon inspection with a flexible inspection scope.

- Ofstead and associates presented the study below at the AORN 2017 Conference<sup>13</sup>. In this study, it was found that sterilized ureteroscopes had high level of visible irregularities, debris, damage and contamination upon inspection with a flexible inspection scope.

**Multisite study on ureteroscope reprocessing effectiveness**  
 Cori L. Ofstead, MSPH<sup>1</sup>; John E. Eiland, RN, MS<sup>2</sup>; Otis L. Heymann, BA<sup>3</sup>; Mariah R. Quick, MPH<sup>1</sup>; Harry P. Wetzler, MD, MSPH<sup>1</sup>  
<sup>1</sup>Ofstead & Associates, Inc., Saint Paul, MN, USA

**Introduction and purpose**

- Contaminated endoscopes, gastroscopes, bronchoscopes, and ureteroscopes have been linked to outbreaks<sup>4,5</sup>
- Design of contaminated endoscopes from area based issues and site base<sup>6,7</sup>
- Functional flaws discovered during reprocessing or reprocessing lead to frequent repairs<sup>8</sup>
- Current guidelines recommend audit and inspection during reprocessing<sup>9</sup>
- Insight sought to answer the following research questions:
  - How much contamination can be observed in a sterile flexible endoscope?
  - How much damage or debris is visible using light inspection?

**Methods**

- Inspection study conducted in two large hospitals
- Site research team
- Audited reprocessing practices
- Obtained samples using surface swabs and a flush broth flush technique
- Performed tests for residual contamination
- Protein, hemoglobin, and albumin (Hematest-SP)<sup>10</sup>
- Microbiological
- Conduct visual inspection of<sup>11</sup>
  - External surface using light inspection and a digital camera
  - Channel and ports using a 3.2 mm fiberoptic bronchoscope

**Results**

- Flexible ureteroscopes (flexible) (n=41):
  - Inspected after 22 cases
  - Reports required after an average of 23 cases for the:
    - Internal leak
    - Insufficient flush quality
    - Broken fibers
    - Blocked channel lumen
- Reprocessing involved:
  - Manual flushing by reprocessing technician
  - Autoclaved with hydrogen peroxide gas
  - Autoclaved with ethylene oxide reprocessing
  - Decontaminated by IR staff
  - Visual inspection by IR and reprocessing staff
  - Doing site re-inspection
  - Contaminated from visible irregularities (Photo 1-3) and contamination (SPH) of endoscopes (Table 1, Figure 1)

Test	Result	Number (%) of endoscopes
Fluorescence	No storage in debris	30 (73%)
Protein	0.4 µg/mL	30 (73%)
Hemoglobin	2.2 µg/mL	1 (8%)
ATP	200 RLU	1 (8%)
Microbial culture	No growth	1 (2%)

**Summary and next steps**

**Sterilized ureteroscopes had high contamination levels, visible damage, and debris**

- Tests conducted to identify flexible endoscopes from:
  - All facilities inspected
  - All facilities contamination above acceptable limits (2 endoscopes)
  - No ATP had positive microbial culture
- Results highlight the need for:
  - Improvement in adherence to guidelines, including:
    - Residue per cleaning by SP-200 to prevent buildup of residue
    - Brochoscope tests that verify cleaning effectiveness
    - Visual inspection with light inspection to identify irregularities
    - More frequent periscope maintenance
    - Reprocessing methods that are proven effective to assess patient safety

**Discussions and acknowledgments**

The study was funded by a grant from the Minnesota State Office of Health Care Regulation and Safety. The authors thank the staff of the participating hospitals for their assistance in conducting the study.

**References**

1. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
2. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
3. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
4. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
5. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
6. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
7. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
8. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
9. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
10. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.
11. Ofstead CL, Eiland JE, Heymann OL, Quick MR, Wetzler HP. Multisite study on ureteroscope reprocessing effectiveness. *J Endourol*. 2017;31(10):1911-1916.

- A poster submitted to IAHCMM for exhibit at the 2018 national conference entitled, “How to get to visually clean: Site and Surface Inspection” detailed the process and requirements for inspecting various types of instruments for residual debris and damage. It recommended that the facility must decide what instruments should be inspected, how often, and with what while keeping in mind current standards and guidelines on the topic. The poster displayed several photos of internal inspection of instruments showing both residual debris and damage that would not qualify the instrument as clean and functional.<sup>14</sup>
- A recent article in HPN by M.A. Drosnack elaborated on the inspection step for instruments<sup>15</sup> and stated that removal of soil and inspection of the instrument for residual debris or damage are critical steps in assessing the quality of reprocessing and whether the item is fit to move on in the process. Items that are deemed to be in disrepair should be taken out of service immediately and those that show signs of having residual soil in or on them either by cleaning verification test or by visual

inspection, should be re-cleaned and re-tested until they pass both the inspection and verification tests.

7. The CMS ASC Infection prevention surveyor worksheet and the CMS Survey and Cert Worksheets state that surveyor should be looking for the following steps to be performed in facilities.
  - a. In section 3.A.6 it is stated, “Flexible endoscopes are inspected for damage and leak tested as part of each reprocessing cycle”.<sup>16</sup>
  - b. Also, in section 3.B.5, the worksheet reads, “Items are thoroughly pre-cleaned according to the manufacturer’s instructions and visually inspected for residual soil prior to sterilization.”<sup>17</sup>

**Example photos demonstrating internal endoscope channels/areas:**




Fig. 1 FIS set up with computer Fig. 2 Internal Channel Fig. 3 Crack in weld at water jet




Fig. 4 Water in Channel




Fig. 5 Debris in Channel

**Example photos of internal suction and shaver lumens:**




Fig. 6 Suction lumens



Fig. 7 Shaver Inspection



Fig. 8 Debris inside shaver




Fig. 9 Shaver debris which tested positive for protein

**Conclusion:**

Therefore, to summarize, the Flexible Inspection Scope (FIS) is a useful tool for inspection of flexible endoscopes, arthroscopy routers, shavers, reamers, flexible endoscopes and other lumened instruments. Use of the FIS borescope provides the end-user a simple verification test to inspect the instrument suction channel, channel openings, distal tip, forceps elevator recess on endoscopes and the internal lumens on other instruments to determine if there are any signs of damage, retained debris and/or moisture remaining in the channel.



## References:

1. ANSI/AAMI ST91:2015. Flexible and semi-rigid endoscope processing in health care facilities.
2. AORN GUIDELINE FOR PROCESSING FLEXIBLE ENDOSCOPES, Revised February 2016 for publication in Guidelines for Perioperative Practice, 2016 edition.
3. ANSI/AAMI ST79:2017 & A1:2010 & A2:2011 & A3:2012 & A4:2013. Comprehensive guide to steam sterilization and sterility assurance in health care facilities.
4. IAHCMM Endoscope Reprocessing Manual, 1<sup>st</sup> Edition. 2017.
5. SGNA POSITION STATEMENT: Management of Endoscopic Accessories, Valves, and Water and Irrigation Bottles in the Gastroenterology Setting. 2018. Available at: <https://www.sgna.org/Portals/0/Management%20Endoscopic%20Accessories%20Valves%20Water%20Irrigation%20bottles.pdf?ver=2018-08-20-141307-367>
6. Arthrex Adaptor Power System TM II (APS II) Shaver Hand pieces DFU-0154r10, Inspection and Maintenance Step 4. [www.arthrex.com](http://www.arthrex.com)
7. STYRKER Shaver Hand Piece – 10000400638 R-2012/10 IFU, Inspection – EN21, Step 9. [www.stryker.com](http://www.stryker.com)
8. FDA Safety Communication. Ongoing Safety Review of Arthroscopic Shavers: FDA Safety Communication. 2015: Archived content. [www.fda.gov](http://www.fda.gov)
9. Ofstead and Associates. 2016. Research poster displayed at AORN 2016. [http://www.ofsteadinsights.com/wp-content/uploads/Ofstead\\_AORN\\_2016\\_poster\\_SUBMITTED.pdf](http://www.ofsteadinsights.com/wp-content/uploads/Ofstead_AORN_2016_poster_SUBMITTED.pdf)
10. Ofstead and Associates, “Assessing Residual Contamination and Damage to Endoscopes over time,” American Journal of Infection Control. 2016
11. Thaker et al., “Inspection of Endoscope Instrument Channels After Reprocessing Using a Prototype Video camera: a Pilot Study,” Abstract presented at DDW National Conference 2017.
12. Ofstead and Associates, “Multisite study on Ureteroscope Reprocessing Effectiveness,” Poster presentation at AORN Surgical Conference. 2017. [http://www.ofsteadinsights.com/wp-content/uploads/Ofstead\\_AORN\\_2017\\_poster.pdf](http://www.ofsteadinsights.com/wp-content/uploads/Ofstead_AORN_2017_poster.pdf)
13. Kovach, S.M., “How to Get to Visually Clean: Site and Surface Inspection,” IAHCMM 2018 National Conference Poster Presentation. <https://www.iahcmm.org/conference-links/poster-gallery/2018-poster-gallery/1305-how-to-get-to-visually-clean.html>
14. Drosnock, M.A., Proper endoscope cleaning requires visual cues — even for the invisible. Oct. 2017. HPN. <https://www.hponline.com/proper-endoscope-cleaning-requires-visual-cues-even-for-the-invisible/>
15. CMS audit worksheet: Exhibit 351 Ambulatory Surgical Center (ASC) INFECTION CONTROL SURVEYOR WORKSHEET (Rev: 142, Issued: 07-17-15, Effective: 07-17-15, Implementation: 07-17-15).
16. Centers for Medicare & Medicaid Services. Hospital Infection Control Worksheet. CMS Survey and Cert Letter 15 12.