

<b>Brand Name of Product</b>	OptiVisor
<b>Generic Name of Product</b>	Headband magnified visor
<b>Product Code Number(s)</b>	26106
<b>Purpose of Product</b>	This product is a precision headband magnifying visor glass lens that improves vision and increases your efficiency. It's used for optical inspection when magnification is needed to visually inspect products to ensure they are properly cleaned.
<b>Range of Applications for Product</b>	Any products that need visual inspection after cleaning.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>Adjustable headband</li> <li>Selection of magnified lens available</li> <li>Lightweight construction</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Wear and use during visual inspection of surgical instruments
<b>Preparation for Use</b>	All of the lens plates are interchangeable. In order to change the lens plate push the fasteners out from the back, inside visor and snap the replacement lens plate.
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>Position the headband by squarely placing on your head with the center of the lens plate at the bridge of the nose.</li> <li>By adjusting the knob at the rear of the unit enough tension can be applied to hold the headband in place.</li> <li>The adjustable pivots on either side of the visor should be tightened with enough tension to hold the lens in front of the eyes as well as up out of the path of vision when not in use.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ol style="list-style-type: none"> <li>The visor should never be raised, lowered, carried, or positioned by means of the lens. Doing so may cause damage both to the lens and the visor.</li> <li>Never put pressure on the lens or area around the lens to position the visor as this may cause the lens to dislocate.</li> <li>This product is not a safety device and should be worn over safety glasses when safety is a concern. The OptiVisor is designed to be worn over prescription eye glasses and not replace them.</li> </ol>
<b>Disposal</b>	

**Reprocessing Instructions**

<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	<ul style="list-style-type: none"> <li>Wash in warm soapy water or lens cleaner.</li> <li>Dry with clean soft cloth (never use abrasive cleaners).</li> </ul>
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	

<b>Storage</b>	
<b>Additional Information</b>	
<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	
<b>Reference Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2018-08-13 msmith