




**Instructions for Use: Pocket Microscope**

<b>Brand Name of Product</b>	Pocket Microscope 20x – 40x
<b>Generic Name of Product</b>	Pocket Microscope
<b>Product Code Number(s)</b>	MG10081-8
<b>Purpose of Product</b>	Perfect carry-with microscope. Ideal for onsite detailed examination of small areas and objects. Useful for quality control inspection to effortlessly identify if the instrument has been properly washed and ready to use.
<b>Range of Applications for Product</b>	Any products that need visual inspection after cleaning.
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Built in light source</li> <li>• Microscope stand allows ambient light to the object under examination.</li> <li>• Camera battery included (3 AG12 batteries)</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	Use during visual inspection of products
<b>Preparation for Use</b>	Remove the Battery Insulate by pulling it outward.
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ul style="list-style-type: none"> <li>• Observing with Specimen Stage               <ol style="list-style-type: none"> <li>1. Place the object of observation onto the Specimen Slide and place a Specimen Slide Cover on top of the object. Insert the Specimen Slide into the Specimen Stage.</li> <li>2. Place the Microscope into the Specimen Stage and press the LED On/Off Button and observe through the Eyepiece.</li> <li>3. Slide the Zoom Lever to adjust the desired degree of magnification (from 20x to 40x)</li> <li>4. Turn the focus wheel gradually until the image is clear and sharp.</li> </ol> </li> <li>• Observing without Specimen Stage               <ol style="list-style-type: none"> <li>1. Place the microscope directly on the top of the object of observation and press the LED On/Off Button.</li> <li>2. Slide the Zoom Lever to adjust the desired degree of magnification (from 20x to 40x)</li> <li>3. Turn the Focusing wheel gradually until the image is clear and sharp.</li> </ol> </li> </ul>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• If the LED light becomes dim, it means that it's time to replace the batteries.               <ol style="list-style-type: none"> <li>1. Unscrew the Battery Door Screw with a screwdriver and remove the Battery Door in the direction shown by the arrow.</li> </ol> </li> </ul>

	<ol style="list-style-type: none"> <li>2. Remove the used-up battery and insert 3 pcs of fresh AG12 button cells into the battery compartment in the right polarity.</li> <li>3. Close the Battery Door and tighten the Battery Door Screw with a screwdriver.</li> </ol>
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	
<b>Reference Documents</b>	Prosys Brochure, Prosys Price List
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2013-04-24 msmith