


Instructions for Use: Pocket Microscope

Brand Name of Product	Pocket Microscope 20x – 40x
Generic Name of Product	Pocket Microscope
Product Code Number(s)	MG10081-8
Purpose of Product	Perfect carry-with microscope. Ideal for onsite detailed examination of small areas and objects. Useful for quality control inspection to effortlessly identify if the item has been properly washed and ready to use.
Range of Applications for Product	Any products that need visual inspection after cleaning.
Key Specifications of Product	<ul style="list-style-type: none"> Built in light source Microscope stand allows ambient light to the object under examination. Camera battery included (3 AG12 batteries)

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product

Description of Use (s)	Use during visual inspection of products
Preparation for Use	Remove the Battery Insulate by pulling it outward.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ul style="list-style-type: none"> Observing with Specimen Stage <ol style="list-style-type: none"> Place the object of observation onto the Specimen Slide and place a Specimen Slide Cover on top of the object. Insert the Specimen Slide into the Specimen Stage. Place the Microscope into the Specimen Stage and press the LED On/Off Button and observe through the Eyepiece. Slide the Zoom Lever to adjust the desired degree of magnification (from 20x to 40x) Turn the focus wheel gradually until the image is clear and sharp. Observing without Specimen Stage <ol style="list-style-type: none"> Place the microscope directly on the top of the object of observation and press the LED On/Off Button. Slide the Zoom Lever to adjust the desired degree of magnification (from 20x to 40x) Turn the Focusing wheel gradually until the image is clear and sharp.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> If the LED light becomes dim, it means that it's time to replace the batteries. <ol style="list-style-type: none"> Unscrew the Battery Door Screw with a screwdriver and remove the Battery Door in the direction shown by the arrow. Remove the used-up battery and insert 3 pcs of fresh AG12 button cells into the

	battery compartment in the right polarity. 3. Close the Battery Door and tighten the Battery Door Screw with a screwdriver.
Disposal	

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instruction	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	
Reference Documents	ProSys™ Brochure, ProSys™ Price List
Customer Service Contact	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2021-03-08 msmith