




**Instructions for Use: DuraSponge™**

<b>Brand Name of Product</b>	DuraSponge™
<b>Generic Name of Product</b>	Instrument Cleaning Sponge
<b>Product Code Number(s)</b>	ESD905201, ESD904101
<b>Purpose of Product</b>	To allow for safe cleaning of sensitive instruments and fixtures, as well as all general surfaces. Tubular sponge is specifically designed for cleaning of flexible and rigid scopes.
<b>Range of Applications for Product</b>	Any instruments that requires immediate removal of soils after procedure.
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Retains 15% more of cleaning solution than standard sponges or paper wipes.</li> <li>• Single-use</li> <li>• Latex Free</li> <li>• Absorbent foam</li> <li>• 2 Styles               <ul style="list-style-type: none"> <li>○ Instrument/Multi-Use-Flat-Dry: 4" W x 5" H</li> <li>○ Tubular/Scope-Dry: 4.5"L x 2.4"D</li> </ul> </li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	<ul style="list-style-type: none"> <li>• Instrument/Multi-Use-Flat-Dry               <ul style="list-style-type: none"> <li>○ Bulked Packaged: 200ct/ 1 case</li> <li>○ Dispenser Box: 200ct/1 case</li> </ul> </li> <li>• Tubular/Scope-Dry               <ul style="list-style-type: none"> <li>○ Bulk Packaged: 100ct/1 case</li> <li>○ Dispenser Box: 100ct/ 1 case</li> </ul> </li> </ul>
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>	
<b>Description of Use (s)</b>	To allow for safe cleaning of sensitive instruments and fixtures, as well as all general surfaces. Tubular sponge is specifically designed for cleaning of flexible and rigid scopes.
<b>Preparation for Use</b>	<ol style="list-style-type: none"> <li>1. Be sure to consult and follow closely the device manufacturer's IFU.</li> <li>2. If utilizing a cleaning agent, be sure to closely follow the detergent manufacturer's IFU for proper preparation of the cleaning solution.</li> </ol>
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ul style="list-style-type: none"> <li>• Moisten sponge with cleaning solution</li> <li>• Wipe down instrument to remove soils and gross debris</li> </ul>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• Be sure to consult the IFU of the device manufacturer to be sure full compliance with their instructions.</li> <li>• Pay special attention to area(s) of the instrument that must be free from moisture.</li> <li>• Pay special attention to chemical compatibility of the medical device and areas of the medical device with the cleaning agent you select.</li> </ul>
<b>Disposal</b>	Treat as biohazard and discard of in compliance with facility guidelines.

**Reprocessing Instructions**

<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Prosys Price Brochure, Prosys Price List
<b>Reference Documents</b>	
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2014-10-01 msmith