



## Instructions for Use: Post Procedure Pre-Treatment Kit

<b>Brand Name of Product</b>	Post Procedure Pre-Treatment Cleaning Kit
<b>Generic Name of Product</b>	Cleaning Kit
<b>Product Code Number(s)</b>	PPPTK-001
<b>Intended Use</b>	For rinsing and wiping of flexible endoscopes at the point of use/bedside.
<b>Range of Applications for Product</b>	Flushing and wiping down of flexible endoscopes after procedures before transporting to the decontamination area.
<b>Key Specifications of Product</b>	<ol style="list-style-type: none"> <li>1. 800 mL of micron filtration water.</li> <li>2. Flat sponge 5”L x 4”W x 0.5”H (127mm x 102mm x 12.5mm), high density, 100% medical grade, biodegradable foam.</li> <li>3. Bag is notched, easy tear open packaging.</li> <li>4. Stable, wide, clear pouch for easy fluid transfer and direct view of the solution.</li> </ol>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	Expires 1 year from the date of manufacture.

Instructions for Using Product	
<b>Description of Use(s)</b>	Flushing and wiping down a flexible endoscope after a procedure following the endoscope manufacturer’s IFU.
<b>Preparation for Use</b>	Open package using the tear notches that are located on both sides of the pouch. Grab the top part of the bag at the notch and tear across to open the bag.
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Wearing appropriate PPE, place a gloved hand inside the bag and squeeze excess fluid out of the flat sponge. Remove the sponge from the pouch.</li> <li>2. Starting at the control section of the scope, place the flat sponge around it and wipe down the length of the insertion tube from the control section to the distal tip to remove any excess bioburden of the surface.</li> <li>3. Take the distal end of the insertion tube, insert it into the pouch and submerge into the solution.</li> <li>4. Suction fluid and air through the endoscope instrument suction channel as instructed by the endoscope manufacturer.</li> <li>5. After several cycles, remove the insertion tube from the solution allowing the solution to drain from the channel by lifting the tube out of the solution so it drains out.</li> <li>6. Follow the endoscope manufacturer’s IFU for any required additional steps for pre-cleaning.</li> <li>7. Place the endoscope into a transport container and transport to the decontamination area.</li> </ol>
<b>Interpretation of Test Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	Store away from direct light. With extended light exposure during storage, the foam sponge may change color.
<b>Disposal</b>	Discard the sponge, bag and fluid using approved methods of disposal.

<b>Reprocessing Instructions</b>	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	
<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	
<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2020-12-01 Suzanne Latta