













**Instructions for Use: Microdenier Wipes**

<b>Brand Name of Product</b>	Ultra Microdenier Wipes
<b>Generic Name of Product</b>	Ultra Microdenier Cleaning Wipes
<b>Product Code Number(s)</b>	TX59
<b>Purpose of Product</b>	For cleaning surfaces requiring low level disinfection.
<b>Range of Applications for Product</b>	The Ultra Microdenier Wipes were specifically developed with an ability to perform over a wide range of applications and settings. Ideal for surfaces that are susceptible to scratching and items that are susceptible to damage by liquid or aerosol cleaners and disinfectants.
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Ultra Microdenier yarn</li> <li>• Capable of removing microbial and other contamination from surfaces without use of liquid solvents</li> <li>• Very tight, stable fabric construction</li> <li>• 100% polyester continuous filament yarn</li> </ul>

<b>Shipping &amp; Storage</b>	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

<b>Instructions for Using Product</b>							
<b>Description of Use (s)</b>	To allow for safe cleaning of sensitive instruments and fixtures, as well as all general surfaces.						
<b>Preparation for Use</b>	1. Be sure to consult and follow closely the device manufacturer's IFU.						
<b>Diagrams (drawings, pictures):</b>	<table border="1"> <tr> <td></td> <td></td> <td></td> </tr> <tr> <td>Fig. 1A unfolded wipe</td> <td>Fig. 1B 1/2 folded wipe</td> <td>Fig. 1C 1/4 folded wipe</td> </tr> </table>  <p>Fig.2</p>				Fig. 1A unfolded wipe	Fig. 1B 1/2 folded wipe	Fig. 1C 1/4 folded wipe
							
Fig. 1A unfolded wipe	Fig. 1B 1/2 folded wipe	Fig. 1C 1/4 folded wipe					
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Follow relevant site protocol (procedures for safety, contamination, etc.) and wear appropriate PPE according to facility policies.</li> <li>2. Fold wipe in mid-air into quarter folds (Fig. 1A–1C). This will produce several clean surface areas and allow better contact with the surface to be wiped.</li> <li>3. When wiping, hold the wiper so that the folded edge is toward the area to be wiped. Hold the unfolded edges in your hand. Group the unfolded edges between thumb and forefinger.</li> <li>4. Wipe in one direction, overlapping wiped area by 10% to 25%.</li> <li>5. Wipe from cleanest to least clean regions of the surface being wiped. Wipe systematically, for example, from top to bottom, far to near. (Fig. 2)</li> <li>6. Keep track of which surfaces have been cleaned and which wiper areas are unused.</li> </ol>						

	7. Always use the cleanest surface of the wiper. If re-wiping use a clean portion of the wiper, not the used wiper area.
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	Be sure to consult the IFU of the device manufacturer to be sure full compliance with their instructions.
<b>Disposal</b>	Treat as biohazard and discard of in compliance with facility guidelines.

<b>Reprocessing Instructions</b>	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	DuraSponge™, Lint Free Wipes
<b>Other Product Support Documents</b>	Prosys Price Brochure, Prosys Price List
<b>Reference Documents</b>	
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2017-08-09 msmith