


## Instructions for Use: Scissor Test Squares

|  |   |
|--|---|
| <b>Brand Name of Product</b>             | Scissor Test Squares  |
| <b>Generic Name of Product</b>           | Scissor Sharpness Test  |
| <b>Product Code Number(s)</b>            | LT-130, RT608, YT606  |
| <b>Purpose of Product</b>                | To test sharpness of scissors   |
| <b>Range of Applications for Product</b> | Scissor sharpness testing   |
| <b>Key Specifications of Product</b>     | <ul style="list-style-type: none"> <li>• LT-130 Latex Free 5 x 5inches– 6.0 gauge</li> <li>• RT608 Latex 216 x 5.5inches– 8.0 gauge</li> <li>• YT606 Latex 216 x 5.5inches – 6.0 gauge</li> </ul> |

| <b>Shipping &amp; Storage</b>                 |  |
|---|--|
| <b>Shipping Conditions &amp; Requirements</b> | N/A  |
| <b>Storage Conditions</b>                     | N/A  |
| <b>Packaging Contents</b>                     | <ul style="list-style-type: none"> <li>• LT-130 quantity: 10 Sheets</li> <li>• RT606 quantity: 1 Roll</li> <li>• YT606 quantity: 1 Roll</li> </ul> |
| <b>Shelf Life</b>                             | N/A  |

| <b>Instructions for Using Product</b>    |  |
|--|--|
| <b>Description of Use(s)</b>             | The test is used to evaluate the sharpness of a scissor.   |
| <b>Preparation for Use</b>               | <ul style="list-style-type: none"> <li>• For total length of the scissor that are less than 5 inches, use the yellow test material or the green LT-130 latex-free option.</li> <li>• For total length of the scissor that are greater than 5 inches, use the red material.</li> </ul>  |
| <b>Diagrams (drawings, pictures)</b>     |   |
| <b>Steps for Use of Product</b>          | <ol style="list-style-type: none"> <li>1. Make at least two cuts in the test material.</li> <li>2. The scissors should cut all the way through the top of the scissor.</li> <li>3. Start cutting at the middle part of the blade cutting clearly through the tip. If you notice any catching or snagging, tag the scissor for repair/sharpening.</li> <li>4. If the scissors are not sharp enough for use, they should be placed in the appropriate repair bin for future inspection by a qualified instrument repair technician according to your medical facility policy.</li> </ol> |
| <b>Interpretation of Results</b>         | If any instrument does not pass the test, it should be given to a qualified instrument repair technician for a closer inspection.  |
| <b>Contraindications of Test Results</b> | N/A  |
| <b>Documentation</b>                     | N/A  |
| <b>Special Warnings and Cautions</b>     | <ul style="list-style-type: none"> <li>• Checking for sharpness in the Central Service Department is part of a quality process.</li> <li>• These are general guidelines to help reduce malfunctioning instruments from being used.</li> <li>• Always remember, if in doubt, place the instrument in the repair bin to have a qualified instrument repair technician check your work.</li> </ul>  |
| <b>Disposal</b>                          | N/A  |

| <b>Reprocessing Instructions</b>       |     |
|--|-----|
| <b>Point of Use</b>                    | N/A |
| <b>Preparation for Decontamination</b> | N/A |
| <b>Disassembly Instructions</b>        | N/A |

|   |     |
|---|-----|
| <b>Cleaning – Manual</b>                    | N/A |
| <b>Cleaning – Automated</b>                 | N/A |
| <b>Disinfection</b>                         | N/A |
| <b>Drying</b>                               | N/A |
| <b>Maintenance, Inspection, and Testing</b> | N/A |
| <b>Reassembly Instructions</b>              | N/A |
| <b>Packaging</b>                            | N/A |
| <b>Sterilization</b>                        | N/A |
| <b>Storage</b>                              | N/A |
| <b>Additional Information</b>               | N/A |

|  |   |
|--|---|
| <b>Related Healthmark Products</b>     | Instrument Repair Tags  |
| <b>Other Product Support Documents</b> | ProSys™ Brochure, ProSys™ Price List  |
| <b>Reference Documents</b>             | N/A   |
| <b>Customer Service Contact</b>        | Healthmark Industries Company, Inc.<br>18600 Malyn Blvd.<br>Fraser, MI 48026<br>1-586-774-7600<br><a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a><br>hmark.com |

2021-05-17 msmith