



**Instructions for Use: Pre-Rinse Shower Unit**

<b>Brand Name of Product</b>	Pre-Rinse Shower Unit
<b>Generic Name of Product</b>	Shower Unit
<b>Product Code Number(s)</b>	300718, 300719, 301189, 301190, 313116, 313323
<b>Purpose of Product</b>	For safe and effective rinsing of gross contaminants from surgical instruments.
<b>Range of Applications for Product</b>	Surgical instruments and other medical devices
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Polished Chrome Plated Brass</li> <li>• Single Hole Deck Mount Pantry Faucet</li> <li>• 1/4 Turn Compression Valves w/ Check Valves</li> <li>• Lever Handles w/ SANIGUARD® Treatment</li> <li>• S/S Flex Hose w/ Santoprene Reinforced Hose</li> <li>• Stainless Steel Strain Relief – Patented</li> <li>• O-Ring Seals - Patented</li> <li>• Standard Full Spray Pattern</li> <li>• 1.2 GPM Spray Face Included</li> <li>• Hold open ring and hose hook included</li> <li>• Commercial Specification Grade</li> <li>• Gaskets and Seals Rated Commercial Hot 180°F</li> <li>• 1/2" NPS Female Inlets</li> <li>• Operating Temperature: 40 - 180°F</li> <li>• Operating Pressure: 15 - 125psi</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	Rinsing of gross contaminates from surgical instruments and equipment.
<b>Preparation for Use</b>	Be sure to consult the IFU of the medical device manufacturer to insure that the device should be exposed to spray washing.
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ul style="list-style-type: none"> <li>• Place instruments soiled and with gross debris under spray unit.</li> <li>• Rinse thoroughly with spray to remove gross debris and soil before further reprocessing.</li> </ul>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• It is usually recommended to use water from the facility's cold water source for the initial rinse in order to avoid denaturing proteins and other organics. Consult the device manufacturer's IFU.</li> </ul>

	<ul style="list-style-type: none"> <li>• The Shower Unit is not a medical product and is therefore not to be used as part of medical treatments.</li> <li>• When in operation, the Shower Unit must not be pointed at parts or orifices of the body because the compressed water or air jet may be harmful.</li> <li>• Be sure to wear proper PPE according to facility and industry guidelines.</li> </ul>
<b>Disposal</b>	

<b>Reprocessing Instructions</b>	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	
<b>Storage:</b>	
<b>Additional Information:</b>	Options include deck mount, splash mount, faucet add-on, wall bracket and an optional 36” long hose.

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Prosys Brochure, Prosys Price List
<b>Reference Documents</b>	Complies with: <ul style="list-style-type: none"> <li>• ANSI A112.18.1M</li> <li>• CSA B125</li> <li>• NSF 61 / 9</li> <li>• ASTM F2324</li> </ul>
<b>Customer Service contact:</b>	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2012-08-13 msmith