




healthmark
INDUSTRIES CO.
health care products
800-521-6224
www.hmark.com

Instructions for Use: Robotic Instrument Guards

Brand Name of Product	Robotic Instrument Guards
Generic Name of Product	Paper Tip Protectors
Product Code Number(s)	PPT-DV-01, PPT-DV-02
Intended Use	Designed to protect the robotic arm tips during sterilization, transport and storage. Prevents the tips of instruments from penetrating the sterile pouch.
Range of Applications for Product	Robotic instruments
Key Specifications of Product	<ul style="list-style-type: none"> • Single use • Two slot sizes available- <ul style="list-style-type: none"> ○ 15mm ○ 30mm • The wide, smooth design allows for easy insertion and removal from slide a sterilization pouch. • Paper protectors are low linting and have excellent wet and dry mechanical properties to withstand the steam sterilization cycle.

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	100/box, paper
Shelf Life	

Instructions for Using Product	
Description of Use(s)	To protect the tip of the robotic instrument during sterilization, transport and storage.
Preparation for Use	
Diagrams (drawings, pictures)	



Steps for Use of Product	<ol style="list-style-type: none"> 1. Fold flap (shorter side) Fig. 1 2. Toward the base (large size) at the score. Fig. 2 3. Tuck the flap into the slotted area. Fig. 3 4. Slide robotic instrument into the first slot, then up through the second slot and back into the pocket formed by the flap and base. Fig. 4 5. Slide robotic arm with paper tip protector into the paper pouch for sterilization. Fig. 5 6. Seal pouch as directed by the instrument pouch IFU.
Interpretation of Test Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	Single-use. Dispose after use in accordance with facility policy.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	<ul style="list-style-type: none"> • Compatible with steam sterilization. • Recommended cycle is 132°C for 4 minutes, or 135°C for 3 minutes in a pre-vac sterilizer. • Paper Material Compatibility <ul style="list-style-type: none"> ○ Steam up to 137°C
Storage	

Additional Information	
Related Healthmark Products	Sterigard Tip Protectors, Vented Instrument Guards, Duo-Vented Guards, Cylindrical Style Instrument Guards
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2020-09-09 Suzanne Latta