



Instructions for Use: Vinyl Instrument Guards

Brand Name of Product	Vinyl Instrument Guards
Generic Name of Product	Tip Protectors
Product Code Number(s)	09-1-011, 09-1-012, 09-1-013, 09-1-014, 09-1-015, 09-1-016, 09-1-017, 09-1-018, 09-1-019, 09-1-020, 09-3-031, 09-3-032, 09-2-033, 09-3-034, 09-3-035, 09-03-037, 09-03-040, 09-07-011, 09-7-012, 09-7-013, 09-7-014, 09-07-15, 09-4-031, 09-4-035, 9-4-036, 09-4-038, 9-4-037, 09-4-032, 09-4-033, 09-4-034
Purpose of Product	Vinyl instrument guards are designed to enhance sterilization and protection for instruments during the sterilization procedure. The vents allow the flow of full-impact sterilization while providing protection to the instruments tips and edges from breaking and/or chipping.
Range of Applications for Product	For trocars, probes and other small tipped instruments.
Key Specifications of Product	<ul style="list-style-type: none"> • Compatible with steam sterilization. • Variety of colors (For Vented and Duo-Vented) • Vented, Duo-Vented and Cylindrical styles

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product

Description of Use (s)	To be placed on small tip surgical instruments
Preparation for Use	Fit each instrument with the appropriate size tip.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place protector tip on appropriate sized instrument prior to sterilization. 2. Package instruments in sterile packaging device (e.g., pouch, wrapped tray, closed container). 3. Place package in sterilizer (steam).
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	

Reprocessing Instructions

Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	<p>They are approved for sterilization via the following methods (Please note, these are minimum cycles:</p> <ul style="list-style-type: none"> • Gravity Steam Autoclave at 121° C for 30 minutes. • Pre-Vac Steam Autoclave at 132 ° C for 3 minutes. • Ethylene Oxide (EtO) • STERRAD® Standard Cycle (approved for web guards only)
Storage:	

Additional Information	
Related Healthmark Products	
Other Product Support Documents	Sterilization Methods For Color Tip Protectors
Reference Documents	Prosys Price List, Prosys Brochure
Customer Service Contact	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2018-01-24 msmith