








Instructions for Use: Window Style Tip Protectors

Brand Name of Product	Window Style Tip Protectors
Generic Name of Product	Instrument Protectors
Product Code Number(s)	15.100.055
Purpose of Product	Designed to protect instrument tips during sterilization, transport and storage.
Range of Applications for Product	For variety of small sizes.
Key specifications of product	<ul style="list-style-type: none"> • Single-use • Paper: <ul style="list-style-type: none"> ○ 250 gram wood free white blanc, allowing excellent performance. • Pouch: <ul style="list-style-type: none"> ○ Transparent pouch made of polyethylene terephthalate (PET) • Ink: <ul style="list-style-type: none"> ○ Non-toxic water-solvent based ink (PH in bet ween 9.0 and 9.5 (at 25°C).

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	Temperature: 10 - 40°C Relative humidity: 20% < RH < 80% Keep away from direct sunlight / keep dry
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	To protect sharp instruments during sterilization, transport and storage.
Preparation for Use	
Diagrams (drawings, pictures):	<div style="display: flex; flex-wrap: wrap;"> <div style="width: 50%; text-align: center;">  <p>1</p> </div> <div style="width: 50%; text-align: center;">  <p>2</p> </div> <div style="width: 50%; text-align: center;">  <p>3</p> </div> <div style="width: 50%; text-align: center;">  <p>4</p> </div> <div style="width: 50%; text-align: center;">  <p>5</p> </div> </div>
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place tip of instrument through top slotted insert. 2. Continue sliding the instrument down the backside of the protector. 3. Tip should emerge through the second slotted insert in front of the protector.

	<p>4. Make sure tip slides inside of window transparent pouch.</p> <p>5. Keep sliding instrument down until its secure.</p>
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	Single use, dispose after use in accordance to laws and facility regulations.

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	<ul style="list-style-type: none"> • Compatible with steam sterilization. Recommended cycle is 4 minutes at 270°F or 3 minutes at 275°F in a pre-vac sterilizer. • Material Compatibility <ul style="list-style-type: none"> ○ Steam up to 145°C ○ Ethylene Oxide
Storage:	
Additional Information:	

Related Healthmark Products	Sterigard Tip Protectors, Vented Instrument Guards, Duo-Vented Guards, Cylindrical Style Instrument Guards
Other Product Support Documents	Prosys Price List, Prosys Brochure
Reference Documents	Paper Point and Needle Protectors comply with the essential requirements of the Medical Devices Directive EC 93/42 - 2007/47/EC - Class 1. They are CE – marked on the product label.
Customer Service contact:	<p>Healthmark Industries Company, Inc</p> <p>33671 Doreka</p> <p>Fraser, MI 48026</p> <p>1-586-774-7600</p> <p>healthmark@hmark.com</p> <p>hmark.com</p>

2017-08-09 msmith