

Brand Name of Product	Humipak®
Generic Name of Product	Transport pouch
Product Code Number(s)	HPSS2035-25, HPSS2035-100, HPSS4243-25, HPSS4243-100, HPSS4267-25, HPSS4267-100, HPSS6577-10, HPSS6577-50
Purpose of Product	To provide a humid environment delaying the drying of contaminants during transport PRIOR to any required cleaning, disinfection, reprocessing and/or sterilization.
Range of Applications for Product	Contaminated items requiring further processing.
Key Specifications of Product	<ul style="list-style-type: none"> • Self-seal strip • A special absorbent layer which when wetted with water creates a humid atmosphere within the pouch • Two layers of transparent film (containing the special absorbent layers) form the pouch • The transparent film, in addition to maintaining the humid atmosphere, allows the contents to be clearly observed • 4 sizes: <ul style="list-style-type: none"> ○ Small, Half Din, Full Din, Double Din <p>Humipak is NOT intended:</p> <ul style="list-style-type: none"> • For medical purposes • To provide any element for cleaning, disinfection or sterilization of reusable medical devices which require cleaning, disinfection or sterilization prior to reuse. • To used as sterile barrier packaging for medical devices.

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	Store in a clean, dry place at room temperature
Packaging Contents	N/A
Shelf Life	5-year shelf life from the date of manufacture

Instructions for Using Product	
Description of Use(s)	Placing items into the Humipak® after use will provide a humid environment delaying the drying of contaminants during transport PRIOR to any required cleaning, disinfection, reprocessing and/or sterilization.
Preparation for Use	Be sure to consult the manufacturer's IFU for pretreatments at the point of use.
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Pick the correct size of Humipak® pouch for the items that are to be kept moist until reprocessing. 2. Add the specified amount of water per the directions found on the Humipak® pouch. 3. Water will wick and be absorbed throughout the Humipak® pouch because of the special absorbent layer between the transparent films. 4. Place items in the pouch. 5. Peel off the protective strip from the adhesive band. 6. Fold over the adhesive strip and press down firmly from the center outwards to ensure a watertight seal. 7. Send to the appropriate area for reprocessing to begin as soon as possible. 8. The Humipak® is NOT a sharps container and is NOT puncture proof. Be sure to follow your facility's policy for handling contaminated reusable sharps in compliance with OSHA Guidelines.
Interpretation of Results	N/A
Contraindications of Test Results	N/A
Documentation	N/A

Special Warnings and Cautions	<ul style="list-style-type: none"> • Items placed in the Humipak[®] should be reprocessed as soon as possible. • Single-use only. • Do not reuse or sterilize the Humipak[®] • Use proper PPE when opening the Humipak[®] • The Humipak[®] is NOT a sharps container and is NOT puncture proof. Be sure to follow your facility's policy for handling contaminated reusable sharps in compliance with OSHA Guidelines.
Disposal	Treat used Humipak [®] as a biohazard and dispose of in accordance with facility's policy for contaminated items.

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	N/A

Related Healthmark Products	SST
Other Product Support Documents	SST Instrument Retrieval Brochure, SST Instrument Retrieval Price List, Humipak [®] Pocket Brochure
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com