



Instructions for Use: Distal Soak

Brand Name of Product	Distal Soak
Generic Name of Product	for transport and storage of laparoscopic & robotic instruments
Product Code Number(s)	50058, 50065, 50069
Intended Use:	The Distal Soak is used during the post-surgical transport and storage of robotic, laparoscopic, or other lumen instruments. The set keeps the distal end moist and facilitates reprocessing of instruments, preventing potential contaminations / soil from drying. Furthermore, the delicate instruments tips are protected during transport.
Range of Applications for Product	Robotic, laparoscopic and other lumen instruments
Key specifications of product	<ul style="list-style-type: none"> • Polypropylene (PP) tubes • Silicone lids • Distal Soak Holder (Sold Separately)

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	<ul style="list-style-type: none"> • The Distal Soak for the transport storage of robotic, laparoscopic or other lumen instruments. • Keeps the distal end moist and enhances the reprocessing of instruments.

Diagrams (drawings, pictures):	
---------------------------------------	--

Steps for Use of Product	<p>Closed Hole Distal Soak</p> <ol style="list-style-type: none"> 1. After the surgical procedure, take the Distal Soak tube out of the packaging unit after. 2. Fill the tube with prepackaged sterile water for irrigation or a cleaning solution that complies with the device manufacturer's IFU. 3. Push the distal end of the instrument through closed hole in the cap. 4. Transport the device to the sterile processing area. 5. When ready to begin further reprocessing, remove the instrument from the Distal Soak. <p>Open Hole Distal Soak</p> <ol style="list-style-type: none"> 1. After the surgical procedure, take the open hole Distal Soak tube out of the packaging unit. 2. Push the distal end of the instrument through the open hole into the <u>empty</u> tube. 3. Inject prepackaged sterile water for irrigation or a cleaning solution that complies with the device manufacturer's IFU. Maximum volume of solution is 15ml. 4. Transport the device to the sterile processing area. 5. When ready to begin further reprocessing, remove the instrument from the Distal Soak.
---------------------------------	---

Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> • Do not use saline solution. • Consult medical device manufacturers IFU for proper precleaning procedures. • Where appropriate PPE at all times when handling contaminated medical devices.
Disposal	Treat used Distal Soak as a biohazard and dispose of in accordance with facility policy for contaminated items.

Reprocessing Instructions	
Point of use:	

Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	
Storage:	
Additional Information:	

Related Healthmark Products	
Other Product Support Documents	
Reference Documents	
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-02-03