

## Instructions for Use: Laser CAAT Tray Tags

<b>Brand Name of Product</b>	Laser CAAT Tray Tags
<b>Generic Name of Product</b>	Tray Tags
<b>Product Code Number(s)</b>	15RBL, 15RGD, 15RGN, 15ROR, 15RPP, 15RRD, 15RSV, 20RBL, 20RGD, 20GRN, 20ROR, 20RPP, 20RRD, 20RSV, 23RCBL, 23RCGD, 23CRD, 23RCSV, 23RCBLK
<b>Purpose of Product</b>	Label trays and other equipment
<b>Range of Applications for Product</b>	
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Anodized aluminum</li> <li>• Sizes: 1.5” round, 2.0” round, 2.5” rectangle</li> <li>• Customized laser engraved</li> <li>• Available in a variety of colors</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use(s)</b>	Label trays and other equipment
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. The Laser CAAT tray tag is attached to the tray or equipment. This should be done using the facility’s specified attachment tool.</li> <li>2. The attachment is threaded through the hole in the tray tag.</li> <li>3. Secure the attachment to the tray or equipment.</li> </ol> <p>Note: Stainless Rings #RIN001 are available for purchase.</p>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• Tray tags should not be in physical contact with the contents of the tray during sterilization cycles.</li> <li>• Do not handle tags until they have cooled to room temperature.</li> </ul>
<b>Disposal</b>	

Reprocessing Instructions	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	Clean with neutral-pH detergent
<b>Cleaning – Automated</b>	<ul style="list-style-type: none"> <li>• Clean with neutral-pH detergent (not alkaline).</li> </ul>
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection, and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	<ul style="list-style-type: none"> <li>• Materially compatible with steam, EO and VH2O2 sterilization.</li> <li>• Recommended cycle time for steam sterilization is 132°C at 4 minutes, 135°C at 3 minutes.</li> </ul>
<b>Storage</b>	
<b>Additional Information</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Labeling Brochure, Labeling Price List
<b>Reference Documents</b>	

<b>Customer Service Contact</b>	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>
---------------------------------	--

2020-08-05 Matt Smith