







Instructions for Use: Transportation Identification Label

| | |
|--|---|
| Brand Name of Product | Transportation Identification Label |
| Generic Name of Product | Tray Tags |
| Product Code Number(s) | AV-52481, AV-52481-RL, AV-52482, AV-52482-RL, HM-CID-613, HM-52483-HLD |
| Intended Use | To identify the status of items during the reprocessing cycle. Label contains pull tabs that indicate if items are clean or dirty. They are single-use. |
| Range of Applications for Product | Case carts or anything used to transport items during reprocessing. |
| Key Specifications of Product | <ul style="list-style-type: none"> • 3.125" x 5.125" perforated tabs • Adhesive backing • Available with and without checklist towards the bottom of the label |

| Shipping & Storage | |
|---|-------------------------------------|
| Shipping Conditions & Requirements | |
| Storage Conditions | Keep in dry, room temperature area. |
| Packaging Conditions | |
| Shelf Life | |

| Instructions for Using Product | | | | | |
|---|--|---|--|-----|-----|
| Description of Use(s) | The Transportation Label can be used on any piece of equipment/cart that is assigned to transport supplies, with the potential to return containing contaminated items. | | | | |
| Preparation for Use | Should be applied to a clean, dry surface. | | | | |
| Diagrams (drawings, pictures) | <table border="1" style="width: 100%; text-align: center;"> <tr> <td></td> <td></td> </tr> <tr> <td>1-3</td> <td>4-5</td> </tr> </table> |  |  | 1-3 | 4-5 |
|  |  | | | | |
| 1-3 | 4-5 | | | | |
| Steps for Use of Product | <ol style="list-style-type: none"> Filling out information on Transportation Identification Label (not all labels have space for information): Upon completion of assembling your transportation cart/bin, fill out the top portion of the label. Using a marker fill out the following (if applicable): <ul style="list-style-type: none"> • Date (When does the procedure take place?) • Case/Procedure (Intended use for the contents) • Room/Facility (Where is the procedure taking place?) • Initials (Who picked/assembled the supplies?) Removing label from the plastic backing: Before placing the label on the transportation cart/bin, slightly bend the bottom ½" of the label on the perforation. This will assist in the removal of the label in decontamination by creating a free piece of the label to grab. Placing the label on the cart/bin: The label should be placed in a visual place for staff to see. This is typically on the door or top of the transportation cart/bin. After the procedure: | | | | |

| | |
|--|--|
| | <p>The staff fills out the transportation checklist on the bottom portion of the label to designate all safety checks have been complete, including staff initials and OR room/facility. This includes marking yes, no or N/A for the following:</p> <ul style="list-style-type: none"> • Linen has been removed • All sharps have been removed • Items are pre-cleaned • Clean and dirty items are separated • Initials and room # <p>NOTE: The staff in decontamination must be able to trace the contaminated items back to a specific location of use or patient in the unfortunate event of a bloodborne pathogen exposure.</p> <p>5. Removing the clean tab on the label:</p> <ul style="list-style-type: none"> • Hands should be covered with a clean pair of gloves for removal of the sticker. • Remove the “Clean” tab by pulling up where it says “Pull Here” in a right to left fashion. • The biohazard symbol and the word “Dirty” will be revealed under the “Clean” tab once it is removed. <p>6. Removing the label in decontamination:</p> <ul style="list-style-type: none"> • Prior to removing the label, clean all the items per your facility policy and empty the transportation cart/bin. • Remove the label by grabbing the bottom area just below the perforation. Pull in a slow motion up and away from the cart. DO NOT PULL FAST. Residue may left behind if pulled too quickly. • Label can be discarded in the garbage. |
| Interpretation of Results | |
| Contraindications of Test Results | |
| Documentation | |
| Special Warnings and Cautions | Do not wash or sterilize with instruments. |
| Disposal | The Transportation Label is a one-time use label and must be discarded after each use in a facility approved waste container according to policy. |

| Reprocessing Instructions | |
|--|--|
| Point of Use | |
| Preparation for Decontamination | |
| Disassembly Instructions | |
| Cleaning – Manual | |
| Cleaning – Automated | |
| Disinfection | |
| Drying | |
| Maintenance, Inspection and Testing | |
| Reassembly Instructions | |
| Packaging | |
| Sterilization | |
| Storage | |
| Additional Information | |

| | |
|--|---|
| Related Healthmark Products | |
| Other Product Support Documents | Labeling Brochure, Labeling Price List |
| Reference Documents | |
| Customer Service Contact | <p>Healthmark Industries Company, Inc. 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com</p> |