









Instructions for Use: Transportation Identification Label

Brand Name of Product	Transportation Identification Label
Generic Name of Product	Tray Tags
Product Code Number(s)	AV52481-250, AV-52481-RL, AV52482-250, AV-52482-RL, HM-CID-613, HM-52483-HLD
Intended Use	To identify the status of items during the transportation cycle. Label contains pull tabs that indicate if items are clean or dirty. They are single-use.
Range of Applications for Product	Case carts, containers, and other transportation methods.
Key Specifications of Product	<ul style="list-style-type: none"> • 3.125" x 5.125" perforated tabs • Adhesive backing • Available with and without checklist towards the bottom of the label

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	Keep in dry, room temperature area.
Packaging Contents	
Shelf Life	

Instructions for Using Product					
Description of Use(s)	The Transportation Label can be used on any piece of equipment/cart that is assigned to transport supplies, with the potential to return containing contaminated items.				
Preparation for Use	Should be applied to a clean, dry surface.				
Diagrams (drawings, pictures)	<table border="1" style="width: 100%; text-align: center;"> <tr> <td></td> <td></td> </tr> <tr> <td>1-3</td> <td>4-5</td> </tr> </table>			1-3	4-5
					
1-3	4-5				
Steps for Use of Product	<ol style="list-style-type: none"> Filling out information on Transportation Identification Label (not all labels have space for information): Upon completion of assembling your transportation cart/bin, fill out the top portion of the label. Using a marker, fill out the following (if applicable): <ul style="list-style-type: none"> • Date (When does the procedure take place?) • Case/Procedure (Intended use for the contents) • Room/Facility (Where is the procedure taking place?) • Initials (Who picked/assembled the supplies?) Removing label from the plastic backing: Before placing the label on the transportation cart/bin, slightly bend the bottom ½" of the label on the perforation. This will assist in the removal of the label in decontamination by creating a free piece of the label to grab. Placing the label on the cart/bin: The label should be placed in a visual place for staff to see. This is typically on the door or top of the transportation cart/bin. After the procedure: 				

	<ul style="list-style-type: none"> • The staff fills out the transportation checklist on the bottom portion of the label to designate all safety checks have been complete. This includes marking yes, no or N/A for the following: <ul style="list-style-type: none"> • Linen has been removed • All sharps have been removed • Items are pre-cleaned • Clean and dirty items are separated • Staff Initials and room # <p>NOTE: In the case of bloodborne pathogen exposure, contaminated items need to be traced back to the specific details of use.</p> <p>5. Removing the clean tab on the label:</p> <ul style="list-style-type: none"> • Hands should be covered with a clean pair of gloves for removal of the sticker. • Remove the ‘Clean’ tab by pulling up where it says ‘Pull Here’ in a right to left fashion. • The biohazard symbol and the word ‘Dirty’ will be revealed under the ‘Clean’ tab once it is removed. <p>6. Removing the label in decontamination:</p> <ul style="list-style-type: none"> • Prior to removing the label, clean all the items per your facility policy and empty the transportation cart/bin. • Remove the label by grabbing the bottom area just below the perforation. Pull in a slow motion up and away from the cart. DO NOT PULL FAST. Residue may be left behind if pulled too quickly. • Label can be discarded in the garbage.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	Do not process with transported items or equipment.
Disposal	The Transportation Label is a single-use label and must be discarded after each use in a facility approved waste container according to policy.

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection and Testing	
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	Labeling Brochure, Labeling Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

