



Instructions for Use: Cynch-Loks

Brand Name of Product	Cynch-Loks
Generic Name of Product	Tamper Evident Seals
Product Code Number(s)	6020, 6323, 2220 LOK
Purpose of Product	To show tampering and/or access to critical equipment and supplies.
Range of Applications for Product	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.
Key specifications of product	<ul style="list-style-type: none"> ● Numbered or Unnumbered ● Cynch-style allows for adjustment to application lengths ● Seal can be broken by hand ● Cannot be resealed once broken ● Variety of colors available ● 2220 LOK designed specifically for lockable round soakers. They communicate “Clean” and “Biohazard” during transportation.

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	100 per box
Shelf Life	

Instructions for Using Product

Description of Use (s)	Seal crash carts, emergency boxes and other areas where evidence of tampering is critical
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<p><u>For Normal Apertures:</u></p> <ol style="list-style-type: none"> 1. Detach bail tip from locking barrel. Pass bail through opening of item to be sealed and through tapered end of locking barrel. 2. Pull bail tightly all the way through. Tug on seal to insure it’s securely locked. <p><u>For Narrow Slot Apertures:</u></p> <ol style="list-style-type: none"> 1. Detach bail tip from locking barrel. Insert tab into narrow slot. 2. Bend tab upwards & pass bail through hole of tab. Twist bail around & insert into tapered end of locking barrel. 3. Pull bail through barrel as far as it will go.
Interpretation of Results	
Contraindications of Test Results	
Documentation	Follow the policy of your facility for recording date, staff and other information critical to management of critical inventory.
Special Warnings and Cautions	If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory.
Disposal	

Reprocessing Instructions

Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	

Reassembly Instructions:	
Packaging:	
Sterilization:	Compatible with steam sterilization. Standard cycle: 270 °/4 minutes
Storage:	
Additional Information:	

Related Healthmark Products	
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-04-05 msmith