


Brand Name of Product	Snap-Loks
Generic Name of Product	Tamper Evident Seals
Product Code Number(s)	7525, 8050
Purpose of Product	To show tampering and/or access to critical equipment and supplies.
Range of Applications for Product	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.
Key Specifications of Product	<ul style="list-style-type: none"> • Numbered and Unnumbered • Seal and broken by hand • Cannot be resealed once broke • Variety of colors available

Shipping & Storage

Shipping Conditions & Requirements	
Storage Conditions	
Packaging Contents	100 per box
Shelf Life	

Instructions for Using Product

Description of Use (s)	Seal crash carts, emergency boxes and other areas where evidence of tampering is critical
Preparation for Use	
Diagrams (drawings, pictures)	 <p>The diagrams illustrate the following steps:</p> <ol style="list-style-type: none"> 1. Inserting the seal into the opening of the object. 2. Pushing the seal through the locking dome. 3. Snapping the seal together to form a tight seal. 4. The final sealed state with the hasp.
Steps for Use of Product	<ol style="list-style-type: none"> 1. Insert seal into closing of two sides of the object being secured. 2. Seal through the locking dome of the seal.Snap item together until Lok is tightly sealed. 3. The seal can be easily opened by breaking by hand by applying pressure at the scored point on the hasp.
Interpretation of Results	
Contraindications of Test Results	

Documentation	Follow the policy of your facility for recording date, staff and other information critical to management of critical inventory.
Special Warnings and Cautions	If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory.
Disposal	

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	
Cleaning – Automated	
Disinfection	
Drying	
Maintenance, Inspection, and Testing	
Reassembly Instructions	
Packaging	
Sterilization	Compatible with steam sterilization. Standard cycle: 132°C (270 °F) for 4 minutes
Storage	
Additional Information	

Related Healthmark Products	
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	
Customer Service Contact	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2021-02-10 msmith