



Instructions for Use: Pull Tight-Loks

Brand Name of Product	Pull Tight-Loks
Generic Name of Product	Tamper Evident Seals
Product Code Number(s)	9393, 5030
Purpose of Product	To show tampering and/or access to critical equipment and supplies.
Range of Applications for Product	Crash carts, emergency boxes, cabinets, and other areas where demonstrating tampering is critical.
Key specifications of product	<ul style="list-style-type: none"> • Numbered or Unnumbered. • Seal and broken by hand. • Design enables sealing items of various dimensions by pulling lock to the optimal seal length. • Cannot be resealed once broken. • Variety of colors available

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	100 per box
Shelf Life	

Instructions for Using Product	
Description of Use (s)	Seal crash carts, emergency boxes and other areas where evidence of tampering is critical
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. To separate seals, hold tab between thumb and forefinger and pull down. 2. Insert the thin end of the Pull Tight seal stem through the locking dome on the side marked "Enter". 3. Pull the stem through the locking dome on the side marked "Enter". Pull it as tight as possible around item being sealed. This will ensure that the seal cannot be broken and re-inserted into the locking mechanism. 4. Utilize the "tug" test by grasping the tab of the seal; gently tug on it to determine closure. If you are able to loosen the seal. IT'S NOT FULLY LOCKED. Remove the seal and replace with a new one.
Interpretation of Results	
Contraindications of Test Results	
Documentation	Follow the policy of your facility for recording date, staff and other information critical to management of critical inventory.
Special Warnings and Cautions	<ul style="list-style-type: none"> • If a seal is broken, assume items have been tampered with or removed. Follow your facility instructions for inspection and replacement of critical inventory. • If the stem is inserted into the wrong side of the tab, the teeth of the locking dome can be damaged and the seal will not work properly. If inserted wrong, remove the seal and begin with a new seal. • Do not re-use a broken seal.
Disposal	

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	

Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	Compatible with steam sterilization. Standard cycle: 270 °/4 minutes
Storage:	
Additional Information:	

Related Healthmark Products	
Other Product Support Documents	Security Products Brochure, Security Products Price List
Reference Documents	
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2015-09-30 msmith