

For more information  
about AccuPoint-HC  
and other cleaning  
verification products,  
call  
800-521-6224 or visit  
HMARK.com

ACCU: 01/2011

## Monitor the Health of Your Facility with



### AccuPoint-HC

### Surface Problems

Hospital acquired infections (HAI) are a significant cause of extended patient stays, rising facility costs and even mortality rates. Independent studies have shown that environmental surfaces are a leading source for HAI's. Verifying that cleaning has been effective has never been easy or quick, until now.



### Rapid Results

Adenosine triphosphate (ATP) is the energy source for all living cells. It is an excellent marker for biological contamination, including key environmental pathogens such as C.dif and S.aureus. The AccuPoint-HC system can analyze a collected sample in 30 seconds and indicate whether a surface is clean.



healthmark  
INDUSTRIES CO.  
health care products

800-521-6224  
www.hmark.com

## Customized for Healthcare

Unlike other systems, AccuPoint-HC has been customized for use in a healthcare setting. This includes the information displayed on the reader as well as the powerful Data Manager software. Track and trend data to spot problem areas, document results, and demonstrate total quality improvement with this powerful software, included with the reader.



## Superior Sampling

The unique design of the collection sampler is superior to traditional swabs. The flat surface comes in complete contact with the area being sampled, insuring a more thorough and accurate result. The self contained cartridge is a complete and easy to use tool whose design minimizes the opportunity for false positive results from poor collection technique.



## AccuPoint-HC is the Answer

Adoption of an auditing program utilizing the AccuPoint-HC removes the subjectivity to evaluating cleaning effectiveness, and will provide actionable data from trending results. Feedback from testing can be used for training of staff, and documentation of facility improvement - always a key focus of Joint Commission, CMS and other auditing bodies.



## Order Today



Code No.	Description	Price
9602HC	AccuPoint HC ATP cleaning validation and tracking system.	\$1,995
9605	AccuPoint Surface Samplers	\$289
	Starter kit: (2 ea 9602HC, 9605 and 1ea 9614)	\$4,595