


















To better clarify differences in Lot numbers used for ChannelCheck items packaged together, please reference the following table:

Product Number	Product Name	Description	Lot numbers (Exterior packaging vs. internal components)	Lot numbers (Test Strips)	Lot number (Water Syringe)
UCC-222	ChannelCheck™	<ul style="list-style-type: none"> <li> 2 Bottles of 50 Test Strips (each)</li> <li> 2 ATS Control Vials</li> <li> 2 Control Color Charts</li> <li> 100 zipper bags</li> <li> 1 IFU</li> </ul>	<p>Exterior packaging Lot number will be different from the internal components.</p> <p>The internal components of the kit have their own Lot number for traceability purposes.</p>	Both bottles of test strips will have the <b>same</b> Lot number.	N/A
UCC-222CP	ChannelCheck™ Convenience Pack (50)	<ul style="list-style-type: none"> <li> 50 Test Strips<sup>Δ</sup></li> <li> 1 ATS Control Vial</li> <li> 50 Zipper bags</li> <li> 50 Pre-filled 10 mL Water Syringes</li> <li> 1 Control Color Chart</li> <li> 1 IFU</li> </ul>	<b>Product Bag*</b> will have the same Lot number as individual test strips within.	Single test strip Lot number <b>matches outside wrap</b> Lot number.	Have separate unique Lot number.
UCC-222TS	ChannelCheck™ Convenience Pack (100)	<ul style="list-style-type: none"> <li> 100 Test Strips<sup>Δ</sup></li> <li> 1 ATS Control Vial</li> <li> 100 Zipper bags</li> <li> 100 Pre-filled 10 mL Water Syringes</li> <li> 1 Control Color Chart</li> <li> 1 IFU</li> </ul>	<b>Product Bag*</b> will have same Lot number as individual test strips within.	Single test strip Lot number <b>matches outside wrap</b> Lot number.	Have separate unique Lot number.

*Δ Individually packaged*

*\*Contains all supplies*

We hope this will give a better understanding to how the Lot numbers are used.

If you have any further questions or concerns, please do not hesitate to [correspond](#).