

HEALTHMARK : ENDOSCOPE SWAB GUIDE

ITEM NUMBER	CHANNEL DIAMETER	SCOPE COLOR
EDH-470, EDP-470	4.7mm-5.0mm	Uncolored
EDH-420, EDP-420	4.0mm-4.6mm	Orange
EDH-350, EDP-350	3.3mm-3.9mm	Orange
EDH-270, EDP-270	2.5mm-3.2mm	Green and Yellow
EDH-200, EDP-200	1.6mm-2.4mm	Purple and Blue
EDH-110, EDP-110	1.0mm-1.5mm	White