

Flexible Endoscopes General Diameter Guide  
for EndoCheck Products

SCOPE TYPE	INSERTION TUBE OUTER DIAMETER	WORKING LENGTH	INSTRUMENT CHANNEL INTERNAL DIAMETER	HEALTHMARK ENDOCHECK PRODUCTS	ENDOCHECK SWAB DIAMETER	PROTIEN DETECTION TEST & SWABS FOR ELEVATOR WIRE	OLYMPUS COLOR GUIDE
<b>ADULTS</b>							
GASTROSCOPE	9.0 mm - 11.4 mm	1030 mm - 1050 mm	2.8 mm - 3.8 mm	EDH-270, EDP-270, EDH-350, EDP-350	2.5-3.2 mm, 3.3-3.9 mm	N/A	
DUODENOSCOPE	10.8 mm - 12.5 mm	1235 mm - 1250 mm	3.2 mm - 4.2 mm	EDH-350, EDP-350, EDH-420, EDP-420	3.3-3.9 mm, 4.0-4.6mm	PT-202 and SWAB-350-12, SWAB-420-12	
COLONOSCOPE	12.9 mm - 13.7 mm	1330 mm - 1680 mm	3.7 mm - 4.2 mm	EDH-350, EDP-350, EDH-420, EDP-420	3.3-3.9 mm, 4.0-4.6mm	N/A	
SIGMOIDOSCOPE	12.8 mm - 13.2 mm	700 mm - 730 mm	3.7 mm - 4.2 mm	EDH-350, EDP-350, EDH-420, EDP-420	3.3-3.9 mm, 4.0-4.6mm	N/A	
ENTEROSCOPE	10.5 mm - 11.7 mm	2200 mm - 2500 mm	2.8 mm - 3.8 mm	EDH-270, EDP-270, EDH-350, EDP-350	2.5-3.2 mm, 3.3-3.9 mm	N/A	
BRONCHOSCOPE	5.7 mm - 6.0 mm	550 mm - 600 mm	2.0 mm - 2.8 mm	EDH-200, EDP-200P, EDH-270, EDP-270	1.6-2.4 mm, 2.5-3.2 mm	N/A	
URETEROSCOPE	2.8 mm - 3.3 mm	670 mm - 700 mm	1.2 mm	EDH-110, EDP-110	1.0-1.5 mm	N/A	White
CYSTOSCOPE	5.4 mm - 5.5 mm	380mm	2.2 mm - 2.4 mm	EDH-200, EDP-200P, EDH-270, EDP-270	1.6-2.4 mm, 2.5-3.2 mm	N/A	
RHINO-LARYNGOSCOPE	2.6 mm - 4.9 mm	300 mm - 365 mm	N/A - 2.0	EDH-200, EDP-200P	1.6-2.4 mm	N/A	
LAPARO-THORACOSCOPE	7.0 mm	270 mm	2.8 mm	EDH-270, EDP-270	2.5-3.2 mm	N/A	
MOBILE AIRWAY SCOPE	4.1 mm - 5.2 mm	600 mm	1.5 mm - 2.6 mm	EDH-110, EDP-110, EDH-200, EDP-200P	1.0-1.5 mm, 1.6-2.4 mm	N/A	White
CHOLEDOCHOSCOPE	2.8 mm - 5.2 mm	380 mm - 700 mm	1.2 mm - 2.2 mm	EDH-110, EDP-110, EDH-200, EDP-200P	1.0-1.5 mm, 1.6-2.4 mm	N/A	White
EUS Linear Scope	11.8 mm - 12.6 mm	1250 mm	2.8 mm - 3.7 mm	EDH-270, EDP-270, EDH-350, EDP-350	2.5-3.2 mm, 3.3-3.9 mm	PT-202, and SWAB-270-12, SWAB-350-12	
EUS Radial Scope	10.5 mm- 11.8 mm	1250 mm	2.2 mm	EDH-200, EDP-200P	1.6-2.4 mm	N/A	
EBUS Scopes	6.2 mm	600 mm	2.2 mm	EDH-200, EDP-200P	1.6-2.4 mm	N/A	
<b>PEDIATRICS</b>							
GASTROSCOPE	5.9 mm - 6.0 mm	1030 mm - 1050 mm	2.0 mm	EDH-200, EDP-200P	1.6-2.4 mm	N/A	
COLONOSCOPE	11.5 mm - 11.6 mm	1680 mm - 1700 mm	3.2 mm - 3.8 mm	EDH-270, EDP-270, EDH-350, EDP-350	2.5-3.2 mm, 3.3-3.9 mm	N/A	
BRONCHOSCOPE	4.4 mm - 5.1 mm	600 mm	2.0 mm	EDH-200, EDP-200P	2.0-2.2 mm	N/A	