

Flexible Endoscopes General Diameter Guide  
for EndoCheck Products

| SCOPE TYPE          | INSERTION TUBE OUTER DIAMETER | WORKING LENGTH    | INSTRUMENT CHANNEL INTERNAL DIAMETER | HEALTHMARK ENDOCHECK PRODUCTS       | ENDOCHECK SWAB DIAMETER | PROTIEN DETECTION TEST & SWABS FOR ELEVATOR WIRE | OLYMPUS COLOR GUIDE |
|---------------------|-------------------------------|-------------------|--------------------------------------|-------------------------------------|-------------------------|--|---------------------|
| <b>ADULTS</b>       |                               |                   |                                      |                                     |                         |  |                     |
| GASTROSCOPE         | 9.0 mm - 11.4 mm              | 1030 mm - 1050 mm | 2.8 mm - 3.8 mm                      | EDH-270, EDP-270, EDH-350, EDP-350  | 2.5-3.2 mm, 3.2-3.9 mm  | N/A  |                     |
| DUODENOSCOPE        | 10.8 mm - 12.5 mm             | 1235 mm - 1250 mm | 3.2 mm - 4.2 mm                      | EDH-420, EDP-420, EDH-350, EDP-350  | 3.3-3.9 mm, 4.0-4.6mm   | PT-202, and SWAB-350-12, SWAB-420-12             |                     |
| COLONOSCOPE         | 12.9 mm - 13.7 mm             | 1330 mm - 1680 mm | 3.7 mm - 4.2 mm                      | EDH-420, EDP-420, EDH-350, EDP-350  | 3.3-3.9 mm, 4.0-4.6mm   | N/A  |                     |
| SIGMOIDOSCOPE       | 12.8 mm - 13.2 mm             | 700 mm - 730 mm   | 3.7 mm - 4.2 mm                      | EDH-420, EDP-420, EDH-350, EDP-350  | 3.3-3.9 mm, 4.0-4.6mm   | N/A  |                     |
| ENTEROSCOPE         | 10.5 mm - 11.7 mm             | 2200 mm - 2500 mm | 2.8 mm - 3.8 mm                      | EDH-270, EDP-270, EDH-350, EDP-350  | 2.5-3.2 mm, 3.2-3.9 mm  | N/A  |                     |
| BRONCHOSCOPE        | 5.7 mm - 6.0 mm               | 550 mm - 600 mm   | 2.0 mm - 2.8 mm                      | EDH-200, EDP-200P, EDH-270, EDP-270 | 1.6-2.4 mm, 2.5-3.2 mm  | N/A  |                     |
| URETEROSCOPE        | 2.8 mm - 3.3 mm               | 670 mm - 700 mm   | 1.2 mm                               | EDH-110, EDP-110                    | 1.0-1.5 mm              | N/A  | White               |
| CYSTOSCOPE          | 5.4 mm - 5.5 mm               | 380mm             | 2.2 mm - 2.4 mm                      | EDH-200, EDP-200P, EDH-270, EDP-270 | 1.6-2.4 mm, 2.5-3.2 mm  | N/A  |                     |
| RHINO-LARYNGOSCOPE  | 2.6 mm - 4.9 mm               | 300 mm - 365 mm   | N/A - 2.0                            | EDH-200, EDP-200P                   | 1.6-2.4 mm              | N/A  |                     |
| LAPARO-THORACOSCOPE | 7.0 mm                        | 270 mm            | 2.8 mm                               | EDH-270, EDP-270                    | 2.5-3.2 mm              | N/A  |                     |
| MOBILE AIRWAY SCOPE | 4.1 mm - 5.2 mm               | 600 mm            | 1.5 mm - 2.6 mm                      | EDH-110, EDP-110, EDH-200, EDP-200P | 1.0-1.5 mm, 1.6-2.4 mm  | N/A  |                     |
| CHOLEDOCHOSCOPE     | 2.8 mm - 5.2 mm               | 380 mm - 700 mm   | 1.2 mm - 2.2 mm                      | EDH-110, EDP-110, EDH-200, EDP-200P | 1.0-1.5 mm, 1.6-2.4 mm  | N/A  |                     |
| EUS Linear Scope    | 11.8 mm - 12.6 mm             | 1250 mm           | 2.8 mm - 3.7 mm                      | EDH-270, EDP-270, EDH-350, EDP-350  | 2.5-3.2 mm, 3.2-3.9 mm  | PT-202, and SWAB-270-12, SWAB-350-12             |                     |
| EUS Radial Scope    | 10.5 mm - 11.8 mm             | 1250 mm           | 2.2 mm                               | EDH-200, EDP-200P                   | 1.6-2.4 mm              | N/A  |                     |
| EBUS Scopes         | 6.2 mm                        | 600 mm            | 2.2 mm                               | EDH-200, EDP-200P                   | 1.6-2.4 mm              | N/A  |                     |
| <b>PEDIATRICS</b>   |                               |                   |                                      |                                     |                         |  |                     |
| GASTROSCOPE         | 5.9 mm - 6.0 mm               | 1030 mm - 1050 mm | 2.0 mm                               | EDH-200, EDP-200P                   | 1.6-2.4 mm              | N/A  |                     |
| COLONOSCOPE         | 11.5 mm - 11.6 mm             | 1680 mm - 1700 mm | 3.2 mm - 3.8 mm                      | EDH-270, EDP-270, EDH-350, EDP-350  | 2.6-3.0 mm, 3.5-3.7 mm  | N/A  |                     |
| BRONCHOSCOPE        | 4.4 mm - 5.1 mm               | 600 mm            | 2.0 mm                               | EDH-200, EDP-200P                   | 2.0-2.2 mm              | N/A  |                     |