



Instructions for Use: Healthmark Watermark™

| | |
|--|--|
| Brand Name of Product | Healthmark Watermark™ |
| Generic Name of Product | Water level guide and temperature gauge |
| Product Code Number(s) | TEMP-WMK |
| Purpose of Product | To indicate the appropriate water volume and temperature in a liquid bath. |
| Range of Applications for Product | When submerging items in a liquid bath |
| Key Specifications of product | |

| Shipping & Storage | |
|---|--|
| Shipping Conditions & Requirements | |
| Storage Conditions | |
| Packaging Contents | |
| Shelf Life | |

| Instructions for Using Product | |
|--|--|
| Description of Use (s) | To gauge the appropriate fill line and temperature of the water. |
| Preparation for Use | <ol style="list-style-type: none"> 1. Be sure the surface for adherence is dry. 2. Read instructions from detergent manufacture and the cleaning equipment or sink. 3. Record the volume of water in the white box on the Healthmark Watermark |
| Diagrams (drawings, pictures) | |
| Steps for Use of Product | <ol style="list-style-type: none"> 1. Determine the appropriate fill line in the tank where the items are to be submerged. 2. Remove adhesive label on the back. 3. Place watermark on water fill line. 4. Add cleaning solution according the instructions of the detergent manufacture. 5. Staff should then be instructed as to the amount of detergent to add to the bath, in accordance with the detergent manufacturer's IFU. |
| Interpretation of Results | |
| Contraindications of Test Results | |
| Documentation | |
| Special Warnings and Cautions | |
| Disposal | |

| Reprocessing Instructions | |
|---|--|
| Point of Use | |
| Preparation for Decontamination | |
| Disassembly Instructions | |
| Cleaning – Manual | The Watermark should be cleaned with the same products used to clean the sink/basin. |
| Cleaning – Automated | |
| Disinfection | The Watermark should be disinfected with the same products as used to disinfect the sink. |
| Drying | |
| Maintenance, Inspection, and Testing | <ul style="list-style-type: none"> • Replace the Healthmark Watermark™ when adhesive back wears off and Watermark begins to peel. • After long term use, the liquid crystals may lose their elasticity and the color change will lose its vibrancy. Change the Watermark should this happen. |
| Reassembly Instructions | |
| Packaging | |
| Sterilization | |
| Storage | |
| Additional Information | |

| | |
|--|--|
| Related Healthmark Products | |
| Other Product Support Documents | ProFormance™ Brochure, ProFormance™ Price List |
| Reference Documents | |

Customer Service Contact

Healthmark Industries Company, Inc.
18600 Malyn Blvd.
Fraser, MI 48026
1-586-774-7600
healthmark@hmark.com
hmark.com