

# ProFormance™ Washer Monitoring Weekly Log Sheet: 4 Level Rack (AWTK-4L)

Date:  Facility:  Washer:  Rack No:  Name:

## Water Quality Aquachek (circle result)

### Target Values

p:

A:

H:




### Comments



### Comments

## Pre-Wash Temperature TempaChek-90

Temperature Should Not Exceed 100°F



Record Result (darken temperatures reached)

### Circle Result

Pass

Fail

### Comments

## Verification of Cleaning Efficacy TOSI

Record Result (circle result)


	Bottom	Middle L	Middle U	Top
0				
1				
2				
3				
4				
5				

### Comments

## Thermal Disinfection TempaChek-170

Target Temp:

Record Result (circle highest temp achieved)



### Comments

Date:	Chamber				Racks						Comments	
	Bottom		Top		Bottom		Middle		Top			
	yes	no	yes	no	yes	no	yes	no	yes	no		
Spray nozzles/arms are free of debris	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Nozzles(holes) properly aligned at target surface (up & down)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
All spray arms are present	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Spray arm spin freely	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Spray arm bushings are intact	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	<input type="checkbox"/>	<input type="checkbox"/>										
Debris screen (in bottom of chamber) is clear of debris	<input type="checkbox"/>	<input type="checkbox"/>										
Instrument rack coupling with manifold properly	<input type="checkbox"/>	<input type="checkbox"/>										
No staining/scaling from detergent, hardwater, etc.	<input type="checkbox"/>	<input type="checkbox"/>										
Detergent/enzyme at sufficient level in container	<input type="checkbox"/>	<input type="checkbox"/>										
Reviewed Cycle Settings (compared to Master)	<input type="checkbox"/>	<input type="checkbox"/>										
Other observations (comment)	<input type="text"/>										<input type="text"/>	
Recommended Actions	<input type="text"/>											