



ProFormanceQA™

Your Central Sterile Quality
Assurance and Infection Control
Documentation Solution

*Reports &
Summaries*

Compiling Summaries...

Choose the desired date range from the 'Summaries' Tab...

That will take you to your Summary Report Menu.

Fill out steps 1 & 2 and click 'Compile'

The screenshot shows a web browser window displaying the ProFormanceQA Summary Report Selection page. The browser's address bar shows the URL: <https://www.blacktiemedical.com/WasherSummaryReport.aspx?dr=4>. The page features the BlackTie logo (Inviting Formality...To The Healthcare Supply Chain) and the healthmark logo (INDUSTRIES CO health care products). The navigation menu includes: Bill-Only Management, Implant Tracking, Contracts & Documents, and ProFormanceQA™. The ProFormanceQA™ menu is expanded, showing a sub-menu with the following options: Setup, Equipment Setup, Conduct Inspections, Results Overview, Reports, and Summaries. The Summaries sub-menu is highlighted with a red box and contains the following date range options: Month to Date, Last Month, Quarter to Date, Last Quarter, Year to Date, and Last Year. The main content area is titled "PerformanceQA Summary Report" and contains three steps: STEP 1: Select Equipment (with a dropdown menu labeled "Select Equipment"), STEP 2: Select Facility, Department and Location (with dropdown menus for "All Facilities", "All Departments", and "All Locations"), and STEP 3: Compile (with a button labeled "Compile" highlighted by a red box). The footer of the page contains the text: "Use of Data - Copyright - Contact Us - Terms of Use - BlackTie PR & Media - Privacy".

The Summary Report...

This is an example of a summary report.

It includes the date range of the report, the locations included, a pass/fail graph, a break down of testing, the dates of all failed tests, and any notes that may be on the log sheets of those failed tests...

ProFormanceQA Summary Report - Windows Internet Explorer

https://www.blacktiemedical.com/WasherSummaryReportPopup.aspx?eq=6&f=SAMPLE FACILITY&d=&l=&dr=0

File Edit View Favorites Tools Help

ProFormanceQA Summary R... ProFormanceQA Summar... x

ProFormanceQA Summary Report - Tunnel Washer

Date Range: 4/1/2009 - 4/16/2009

Results include the following:

Facility	Department	Location
SAMPLE FACILITY	SAMPLE DEPT.	SAMPLE LOCATION

Key Quality Indicator: TOSI

Total # of KQI Scheduled:	32
Total # of KPI Run:	4
% Run vs Scheduled:	12%
# Pass:	3
# Fail:	1
% Pass vs Fail:	75%

Dates of Failed Tests:

Date	Facility	Department	Location	Equipment ID
4/16/2009	SAMPLE FACILITY	SAMPLE DEPT.	SAMPLE LOCATION	Tunnel Washer 456

Inspector Notes:

Pass/Fail

A 3D pie chart titled "Pass/Fail" showing the distribution of test results. The chart is divided into two segments: a large green segment representing "Pass - 75%" and a smaller red segment representing "Fail - 25%".

Pass - 75%
Fail - 25%

Print Close

Done Internet 100%

Generating Reports...

Choose the desired equipment from the "Reports/Chart Results" Tab...
That will take you to your Report Menu.

Choose the date range of your report.

Choose the facilities, departments, locations, and inspectors that you would like to include.

Click "Generate"

The screenshot shows a web browser window displaying the BlackTie ProFormanceQA™ application. The browser title is "Tunnel Washer Report - Windows Internet Explorer" and the address bar shows "https://www.blacktiemedical.com/TunnelWasherReport.aspx". The application header includes the BlackTie logo, navigation links (Home, About Us, Contact, Log Out), and the user name "User: SAMPLE SAMPLE". The main navigation menu includes "Bill-Only Management", "Implant Tracking", "Contracts & Documents", and "ProFormanceQA™". The "ProFormanceQA™" menu is expanded, showing a sub-menu with "Reports" and "Summaries". The "Reports" sub-menu is further expanded, showing a list of report types: Hospital Staff, Auto Endoscope Reprocessor, Cart Report, Channel Check Report, EndoCheck Report, HemoCheck Report, Lumen Report, Multi Rack Washer Report, Spot Check Report, Tunnel Washer Report, Ultrasonic Report, Washer Staff, and Water Quality. The "Tunnel Washer Report" option is highlighted. Below the navigation menu, the "Tunnel Washer Report" form is visible, with fields for "Date" (Quick, From: 04/09/2009, To: 04/16/2009), "Facility" (All), "Department" (All), "Location" (All), and "Inspector" (All). A "Generate" button is highlighted with a red box.

The Report...

This is an example of a report.

It includes graphs for water temperature and a line detail report below.

The line detail report shows all data entered on the log sheets.

Tunnel Washer Report - Windows Internet Explorer

https://www.blacktiemedical.com/TunnelWasherReport.aspx

File Edit View Favorites Tools Help

Tunnel Washer Report

Bill-Only Management Implant Tracking Contracts & Documents **ProFormanceQA™**

Setup Equipment Setup Conduct Inspections Results Overview Reports Summaries

Tunnel Washer Report

Quick: From: To:

Date: This Month Or 04/09/2009 04/16/2009

Facility: SAMPLE FACILITY Department: SAMPLE DEPT. Location: SAMPLE LOCATION

Inspector: All Generate

TempChek-90 (F)

Zoom

TempChek-170 (F)

Zoom

Date	First Name	Last Name	Location	TempaChek-90 (F)	TempaChek-170 (F)	SonoCheck	SonoCheck Placement	TOSIDegassed?	Nozzles Free Pre Bottom	Nozzles Debris Pre Top	Nozzles Free Sides	Nozzles Debris Wash Bottom	Nozzles Free Wash Top	Nozzles Debris Sides	Nozzles Free Bottom	Nozzles Debris Bottom	Arms Present Bottom	Arms Present Pre Top	Arms Present Sides	Arms Present Wash Bottom	Arms Present Wash top	Arms Present Disinfect Sides	Arms Present Disinfect Bottom	Arms Present Top	Arms Present Sides	Arms Present Pre
4/16/2009	JOHN	DOE	SAMPLE LOCATION	100	180	Yellow	Center	0	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
4/16/2009	JOHN	DOE	SAMPLE LOCATION	100	190	Green	Center	3	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
4/16/2009	JOHN	DOE	SAMPLE LOCATION	100	180	Yellow	Center	0	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
4/16/2009	JOHN	DOE	SAMPLE LOCATION	100	180	Yellow	Center	0	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes

Use of Data - Copyright - Contact Us - Terms of Use

Congratulations!

You have completed
ProFormanceQA
Reports & Summaries Training.

If you have any questions or comments,
please contact Alex Kirk at Healthmark Industries.

Phone: (800)521-6224

Email: alex@hmark.com

Thank You!