

# ProFormance™ Log Sheet: 3 Level Rack (WTK-3L)

Date:  Facility:  Washer:  Rack No:  Name:

	Chamber				Racks				Comments
	Bottom		Top		Bottom		Top		
	yes	no	yes	no	yes	no	yes	no	
Spray nozzles/arms are free of debris	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Nozzles(holes) properly aligned at target surface (up & down)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
All spray arms are present	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Spray arm spin freely	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	yes	no							
Debris screen (in bottom of chamber) is clear of debris	<input type="checkbox"/>	<input type="checkbox"/>							
Instrument rack coupling with manifold properly	<input type="checkbox"/>	<input type="checkbox"/>							
No staining/scaling from detergent, hardwater, etc.	<input type="checkbox"/>	<input type="checkbox"/>							
Detergent/enzyme at sufficient level in container	<input type="checkbox"/>	<input type="checkbox"/>							
Reviewed Cycle Settings (Compared to Master)	<input type="checkbox"/>	<input type="checkbox"/>							

				Comments
<b>TOSI</b>	Bottom	Middle	Top	<div style="background-color: #ff9933; height: 20px; margin-bottom: 5px;"></div> <div style="background-color: #ff9933; height: 20px; margin-bottom: 5px;"></div> <div style="background-color: #ff9933; height: 20px; margin-bottom: 5px;"></div> <div style="background-color: #ff9933; height: 20px; margin-bottom: 5px;"></div> <div style="background-color: #ff9933; height: 20px; margin-bottom: 5px;"></div> <div style="background-color: #ff9933; height: 20px;"></div>
Record Result (circle result)	0	1	2	
	3	4	5	
				
				
				