



ProFormance™ Log Sheet: 4 Level Rack (WTK-4L)

Date: Facility: Washer: Rack No: Name:

	Chamber				Racks						Comments
	Bottom		Top		Bottom		Middle		Top		
	yes	no	yes	no	yes	no	yes	no	yes	no	
Spray nozzles/arms are free of debris	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
Nozzles(holes) properly aligned at target surface (up & down)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>
All spray arms are present	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>
Spray arm spin freely	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>
	yes	no	Comments								
Debris screen (in bottom of chamber) is clear of debris	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>								
Instrument rack coupling with manifold properly	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>								
No staining/scaling from detergent, hardwater, etc.	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>								
Detergent/enzyme at sufficient level in container	<input type="checkbox"/>	<input type="checkbox"/>	<input style="width: 100%; height: 20px;" type="text"/>								

		Bottom	Middle1	Middle2	Top	Comments
TOSI						
Record Result (circle result)	0					<input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/> <input style="width: 100%; height: 20px;" type="text"/>
	1					
	2					
	3					
	4					
	5					