


## Instructions for Use: LTA Disposable External Brush-Rigid Fibers

<b>Brand Name of Product</b>	LTA Disposable External Brush
<b>Generic Name of Product</b>	Disposable Brush
<b>Product Code Number(s)</b>	SU-510, SU-520
<b>Intended Use</b>	Designed for preliminary cleaning of hard surface items prior to further processing.
<b>Range of Applications for Product</b>	Products that require precision cleaning of devices with hard external surfaces.
<b>Key Specifications of Product</b>	<ul style="list-style-type: none"> <li>• Single-use brushes are designed with latex-free thermoplastic handles and nylon bristles.</li> <li>• Medium fiber style brush is colored yellow with white bristles and is recommended for routine cleaning.</li> <li>• Rigid fiber style brush is colored green with white bristles and is recommended for intensive cleaning.</li> </ul>

Shipping & Storage	
<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Contents</b>	
<b>Shelf Life</b>	

Instructions for Using Product	
<b>Description of Use(s)</b>	Designed for preliminary cleaning of hard surface items.
<b>Preparation for Use</b>	Before each use, check the integrity and the condition of the brush. Replace any worn or damaged brush.
<b>Diagrams (drawings, pictures)</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>1. Consult the manufacturer's IFU for the item you are cleaning and follow those instruction closely. Pay close attention to recommendations on the material composition of the brush trim.</li> <li>2. Visually inspect box lock area and serrations for any debris / residual / bioburden / blood / protein.</li> </ol>
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	<ul style="list-style-type: none"> <li>• External brushes are intended for use as cleaning tools only. They are not for use in patient care or aseptic techniques.</li> <li>• Wear proper PPE in compliance with facility policy.</li> </ul>
<b>Disposal</b>	Because these brushes are used to clean patient used items, it is recommended to dispose of them in a biohazard container.

Reprocessing Instructions	
<b>Point of Use</b>	
<b>Preparation for Decontamination</b>	
<b>Disassembly Instructions</b>	
<b>Cleaning – Manual</b>	

<b>Cleaning – Automated</b>	
<b>Disinfection</b>	
<b>Drying</b>	
<b>Maintenance, Inspection and Testing</b>	
<b>Reassembly Instructions</b>	
<b>Packaging</b>	
<b>Sterilization</b>	
<b>Storage</b>	
<b>Additional Information</b>	Intended for single-use. A “use” should be defined by facility policy.
<b>Related Healthmark Products</b>	Brush Cleaning Tools
<b>Other Product Support Documents</b>	ProSys™ Brochure, ProSys™ Price List
<b>Reference Documents</b>	
<b>Customer Service Contact</b>	Healthmark Industries Company, Inc 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 <a href="mailto:healthmark@hmark.com">healthmark@hmark.com</a> <a href="http://hmark.com">hmark.com</a>

2021-01-27 Suzanne Latta