



**Instructions for Use: Clamps & Spreaders**

<b>Brand Name of Product</b>	Clamps & Spreaders
<b>Generic Name of Product</b>	Clamps & Spreaders
<b>Product Code Number(s)</b>	20411, 20408, 20406, 20412, 20407
<b>Purpose of Product</b>	To keep surgical instruments in an open position while being cleaned.
<b>Range of Applications for Product</b>	<ul style="list-style-type: none"> <li>• Clamps ensure that joints of jaws of the bone forceps remain open.</li> <li>• Spreaders are for instruments with ring handles that cannot be dismantled. They can be cleaned more efficiently in the washer-disinfector. The jaws remain open during the cleaning process, thus assuring thorough cleaning.</li> </ul>
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>• Stainless Steel</li> <li>• Silicone Covering</li> </ul>

**Shipping & Storage**

<b>Shipping Conditions &amp; Requirements</b>	
<b>Storage Conditions</b>	
<b>Packaging Conditions</b>	
<b>Shelf Life</b>	

**Instructions for Using Product**

<b>Description of Use (s)</b>	When instruments are closed, they are not cleaned properly. Water and detergent cannot reach all of the surfaces of the instrument. This leads to an increase chance of contamination.
<b>Preparation for Use</b>	Make sure you select the appropriate size clamp before using on instrument.
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ul style="list-style-type: none"> <li>• Open surgical instrument to the appropriate width for cleaning.</li> <li>• Insert ends of clamps or spreaders through instrument loops. This will help keep medical device in an open position for cleaning.</li> <li>• Continue normal cleaning procedure according to instrument manufacture guidelines.</li> </ul>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	AAMI, AORN, device manufactures and others strongly recommend that instruments be kept in an open position when cleaned, to insure maximum exposure to the cleaning solution.
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

**Reprocessing Instructions**

<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	
<b>Cleaning – Automated:</b>	Can be cleaned through the automatic washer with surgical instruments.
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	Materially compatible for Steam and EO sterilization.
<b>Storage:</b>	
<b>Additional Information:</b>	

<b>Related Healthmark Products</b>	
<b>Other Product Support Documents</b>	Prosys Brochure, Prosys Price List
<b>Reference Documents</b>	

**Customer Service contact:**

Healthmark Industries Company, Inc  
33671 Doreka  
Fraser, MI 48026  
1-586-774-7600  
[healthmark@hmark.com](mailto:healthmark@hmark.com)  
[hmark.com](http://hmark.com)

2017-08-09 msmith