

Instructions for Use: Silicone Heat Glove

Brand Name of Product	Silicone Heat Glove
Generic Name of Product	Heat Glove
Product Code Number(s)	09484, 09368
Intended Use	For the safe handling of hot trays and racks from autoclaves
Range of Applications for Product	Tray, racks and items which are hot that are being removed from the autoclave
Key Specifications of Product	<ul style="list-style-type: none"> Made of Silicone Protects from surface temperatures ranging from -40°F to +464°F

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	Stored at room temperature. Do not expose to direct sunlight.
Packaging Contents	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	
Preparation for Use	
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> Insert hands completely into the glove. Be sure that the glove is completely covering your hand for maximum protection.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	<ul style="list-style-type: none"> The gloves can stand up to temperatures of 240°C (464°F). Gloves do not protect against chemicals.
Disposal	

Reprocessing Instructions	
Point of Use	
Preparation for Decontamination	
Disassembly Instructions	
Cleaning – Manual	<ul style="list-style-type: none"> Detergent (pH < 10). Consult the detergents' IFU for the recommended concentration, exposure times and temperature. It is recommended to use demineralized water for the final rinse.
Cleaning – Automated	<ul style="list-style-type: none"> Acidic / neutral / alkaline < pH 10. Consult the detergents' IFU for the recommended concentration, exposure times and temperature. Consult the IFU of the washer manufacturer for correct cycle settings. Do not clean in an ultrasonic cleaner.
Disinfection	
Drying	
Maintenance, Inspection, and Testing	<ul style="list-style-type: none"> Visually inspect gloves for usability and cleanliness. If gloves are damaged, they no longer provide protection and must be replaced.
Reassembly Instructions	
Packaging	
Sterilization	
Storage	
Additional Information	

Related Healthmark Products	Seamless Autoclave Gloves, Terry Cloth Autoclave Gloves, Shark Tooth Heat Mitts
Other Product Support Documents	ProSys™ Brochure, ProSys™ Price List
Reference Documents	

Customer Service Contact

Healthmark Industries Company, Inc
18600 Malyn Blvd.
Fraser, MI 48026
1-586-774-7600
healthmark@hmark.com
hmark.com

2019-07-03 msmith

Rev. B