

System Code Number	Description	Glutaraldehyde	Ortho-phthalaldehyde	Peracetic Acid	Enzymatic Detergent	Neutral Detergent	Hydrogen Peroxide Wipes	Liquid Hydrogen Peroxide	Alcohol Wipes	Bleach Wipes	Additional Comments
FIS-001	Flexible Inspection Scope Camera						Y		Y	Y	
FIS-003	Flexible Inspection Scope Camera						Y		Y	Y	
FIS-004	2.4mm Video Borescope					Y	Y		Y		Isopropyl Alcohol
FIS-005, FIS-005SK	Flexible Inspection Scope Kit	Y	Y	Y	Y	Y			Y		Camera Cable- Isopropyl Alcohol, Can Not Be Submerged
MICR007	USB Microscope								Y		Ethyl Alcohol
FBR 001, FBR 002	Fiberscope and LED Light Source Adap.	Y	Y		Y	Y			Y		LED Light Source Adapter- Enzymatic and Neutral pH solutions