

	FIS-S001	FIS-003	FIS-004	FIS-005, FIS-005SK	MICR007	FBR 001, FBR 002
Description	Flexible Inspection Scope Camera	Flexible Inspection Scope Camera	2.4mm Video Borescope	Flexible Inspection Scope Kit	USB Microscope	Fiberscope and LED Light Source Adap.
Glutaraldehyde				Y		Y
Ortho-phthalaldehyde				Y		Y
Peracetic Acid				Y		
Quaternary Ammonia Wipes	Y	Y	Y	Y	Y	Y
Enzymatic Detergent				Y		Y
Neutral Ph Detergents			Y	Y		Y
Hydrogen Peroxide Wipes	Y	Y	Y			
Liquid Hydrogen Peroxide				Y		Y
Alcohol Wipes	Y	Y	Y	Y	Y	Y
Bleach Wipes	Y	Y				
Sodium Borate				Y		Y
Sodium Xylene Sulfonate				Y		Y
Tetrapotassium Pyrophosphate				Y		Y
Trisodium Nitrilotriacetate				Y		Y
Ethylene oxide				Y		Y
STERRAD® Systems (Standard Cycle)				Y		Y
Steris System 1E				Y		Y
STERIS V-PRO® Low Temperature Sterilization Systems (NonLumen Cycle)				Y		Y
Additional Comments			Isopropyl Alcohol	Camera Cable- Isopropyl Alcohol, Can Not Be Submerged	Ethyl Alcohol	LED Light Source Adapter- Enzymatic and Neutral pH solutions