



Instructions for Use: Multi-Magnifier

Brand Name of Product	Multi-Magnifier
Generic Name of Product	Multi-Magnification Magnifier
Product Code Number(s)	MAG-310
Purpose of Product	Small size that makes it convenient to carry at all times. Designed to allow for easy and thorough inspection of surgical instruments. All in one magnifier magnifies by 3x, 10x and 55x
Range of Applications for Product	Surgical instruments that need up close visual inspection for possible damage or residual soiling after reprocessing.
Key specifications of product	<ul style="list-style-type: none"> • 3x, 10x and 55x magnification • Lightweight and durable design • Built in lighting • Includes battery

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	Close-up visual inspection of instruments after reprocessing
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ol style="list-style-type: none"> 1. After reprocessing; use pocket magnifier at optical inspection stations. 2. Select lens for use during optical inspection by amount needed. Magnifies by 3x, 10x and 55x. 3. Press button on side to turn on light source. Each time you press the button the light source will change light intensity and location of the lights that illuminate.
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	In order to prevent scratching, clean the lens using lens cleaner and a soft cloth.
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	
Storage:	

Additional Information:	
Related Healthmark Products	
Other Product Support Documents	
Reference Documents	ProSys Brochure, ProSys Price List
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2016-12-19 Ashlyne Basile