



healthmark
INDUSTRIES CO.
health care products
800-521-6224
www.hmark.com

Instructions for Use: Wire Mesh & Larger Washer Baskets

Brand Name of Product	Wire Mesh and Larger Washer Baskets
Generic Name of Product	Washer baskets for medium and larger items
Product Code Number(s)	75.0005.3, 75.0004.6, 75.75.0005.9, 75.0005.8, 75.0004.8, 75.0004.9, 75.0091.2, 75.0091.1, 72.1400.5, 72.1400.6
Purpose of Product	To hold parts and instruments when being washed inside of an automated instrument washer, preventing damage and loss.
Range of Applications for Product	Instruments of from small to larger sizes.
Key specifications of product	<ol style="list-style-type: none"> 1. Manufactured from stainless steel. 2. Lid for safe containment of small parts. (Available on the Mesh Washer Basket only) 3. Designed to allow thorough exposure to cleaning and rinsing solution within an automated instrument washer.

Shipping & Storage	
Shipping Conditions & Requirements	
Storage Conditions	
Packaging Conditions	
Shelf Life	

Instructions for Using Product	
Description of Use (s)	Place in the washer basket small surgical instruments and parts that often get lost or damaged inside of an automated instrument washer.
Preparation for Use	
Diagrams (drawings, pictures):	
Steps for Use of Product	<ul style="list-style-type: none"> • Place instruments in basket • Close lid of basket (If applicable) • Place into instrument washer for cleaning
Interpretation of Results	
Contraindications of Test Results	
Documentation	
Special Warnings and Cautions	
Disposal	

Reprocessing Instructions	
Point of use:	
Preparation for decontamination:	
Disassembly Instructions:	
Cleaning – Manual:	
Cleaning – Automated:	
Disinfection:	
Drying:	
Maintenance, inspection, and testing:	
Reassembly Instructions:	
Packaging:	
Sterilization:	Compatible with steam sterilization. Recommended cycle is 4 minutes at 270°F or 3

	minutes at 275°F in a pre-vac sterilizer.
Storage:	
Additional Information:	

Related Healthmark Products	Nic-Nac Washer Baskets and Balls
Other Product Support Documents	Prosys Brochure, Prosys Price List
Reference Documents	
Customer Service contact:	Healthmark Industries Company, Inc 33671 Doreka Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com

2014-07-08 msmith