


| | |
|--|---|
| Brand Name of Product | Face Shield with Drape |
| Generic Name of Product | Face Shield with Drape |
| Product Code Number(s) | SHL-002 |
| Intended Use | Designed to help protect the wearer from splash. |
| Range of Applications for Product | Areas where protection against splash is indicated. |
| Key Specifications of Product | <ul style="list-style-type: none"> • Latex-free elastic headband • Antistatic face shield • Face shield dimensions: 13 in. (H) x 9 in. (W) • 7 mil thick clear polyester (PET) • One inch (25.4 mm) thick-polyurethane brow foam • Drape dimensions- 12 inches in Length (300 mm) x 12 inches in Width (300 mm) • Independently tested and meets the minimum liquid barrier performance criteria as specified in ASTM F2100. |

| | |
|---|---|
| Shipping & Storage | |
| Shipping Conditions & Requirements | N/A |
| Storage Conditions | Product must be stored out of sunlight and in a clean, dry (between 30% rh and 60% rh area maintained at a temperature between 10 °C and 30 °C (50 °F and 86 °F). |
| Packaging Contents | N/A |
| Shelf Life | N/A |

| | |
|--|--|
| Instructions for Using Product | |
| Description of Use(s) | Designed to help protect the wearer from splash. |
| Preparation for Use | N/A |
| Diagrams (drawings, pictures) |  <p>Areas of Shield in Direct Contact with Wearer:</p> <ul style="list-style-type: none"> • Forehead • Above Ears • Around Back of Head • Chest, when drape is tucked |
| Steps for Use of Product | <ol style="list-style-type: none"> 1. Place foam front against your forehead. 2. Secure elastic band around the back of your head. 3. Tuck the drape into torso garment so that it covers the front facing neck area. |
| Interpretation of Test Results | N/A |
| Contraindications of Test Results | N/A |
| Documentation | N/A |
| Special Warnings and Cautions | <ul style="list-style-type: none"> • This face shield or any PPE is to be used according to your specific facility policy. |

| | |
|-----------------|--|
| | <ul style="list-style-type: none"> • The face shield alone will not prevent infection from microbes or viruses nor is it effective against radiation. • Healthcare facilities are responsible for decisions regarding appropriate PPE for their personnel. |
| Disposal | Single-use |

| Reprocessing Instructions | |
|---|---|
| Point of Use | N/A |
| Preparation for Decontamination | N/A |
| Disassembly Instructions | N/A |
| Cleaning – Manual | N/A |
| Cleaning – Automated | N/A |
| Disinfection | N/A |
| Drying | N/A |
| Maintenance, Inspection, and Testing | N/A |
| Reassembly Instructions | N/A |
| Packaging | N/A |
| Sterilization | N/A |
| Storage | Product must be stored out of sunlight and in a clean, dry (between 30% rh and 60% rh) area maintained at a temperature between 10 °C and 30 °C (50 °F and 86 °F). |
| Additional Information | Independent test results are available upon request. |
| Related Healthmark Products | N/A |
| Other Product Support Documents | SST Brochure, SST Price List |
| Reference Documents | N/A |
| Customer Service Contact | Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com |