


Instructions for Use: Premium Face Shield

Brand Name of Product	Premium Face Shield
Generic Name of Product	Face Shield
Product Code Number(s)	SHL-003
Intended Use	Designed to protect the wearer from splash.
Range of Applications for Product	Areas where protection against splash is indicated.
Key Specifications of Product	<p>Shield</p> <ul style="list-style-type: none"> • Anti-fog PET Plastic • Shield- 19.5” wide x 8.5” tall • Clear plastic- 0.020” x 20” x 10” <p>Foam</p> <ul style="list-style-type: none"> • 1” thick ultra-comfort polyurethane • Hypoallergenic foam band can absorb sweat and provide enough room for eyeglasses or safety goggles <p>Head Band</p> <ul style="list-style-type: none"> • 1 ¼” black knitted elastic band • 35% elastic • 65% polyester <p>Has been independently tested and meets the liquid barrier criteria as specified in ASTM F2100.</p>

Shipping & Storage	
Shipping Conditions & Requirements	N/A
Storage Conditions	N/A
Packaging Contents	N/A
Shelf Life	N/A

Instructions for Using Product	
Description of Use(s)	Designed to protect the wearer from splash.
Preparation for Use	N/A
Diagrams (drawings, pictures)	
Steps for Use of Product	<ol style="list-style-type: none"> 1. Place foam front against your forehead. 2. Secure elastic band around the back of your head.
Interpretation of Test Results	N/A
Contraindications of Test Results	N/A

Documentation	N/A
Special Warnings and Cautions	<ul style="list-style-type: none"> • The face shield is not a barrier to microbes or viruses. • The use of this or any PPE is to be used according to your specific facility policy. • Healthcare facilities are responsible for decisions regarding appropriate personal protective equipment for their personnel.
Disposal	Single-use

Reprocessing Instructions	
Point of Use	N/A
Preparation for Decontamination	N/A
Disassembly Instructions	N/A
Cleaning – Manual	N/A
Cleaning – Automated	N/A
Disinfection	N/A
Drying	N/A
Maintenance, Inspection, and Testing	N/A
Reassembly Instructions	N/A
Packaging	N/A
Sterilization	N/A
Storage	N/A
Additional Information	Independent test results available upon request
Related Healthmark Products	N/A
Other Product Support Documents	PPE Brochure, PPE Price List
Reference Documents	N/A
Customer Service Contact	Healthmark Industries Company, Inc. 18600 Malyn Blvd. Fraser, MI 48026 1-586-774-7600 healthmark@hmark.com hmark.com